# LED ILLUMINATION

## OPERATING INSTRUCTION


ryf ag



Ryf AG  
Bettlachstrasse 2  
2540 Grenchen  
tel 032 654 21 00  
fax 032 654 21 09


[www.ryfag.ch](http://www.ryfag.ch)


(A)




(B)

(C)

(D)

(E)

- (A) 5 Plug for 4
- 7 Fixing screw
- (B) 5 Plug for 4
- 6 Gooseneck
- (C) 8 Plug for mains cable
- 9 Plug for 3
- (D) 1 On-Off Switch
- 2 Brightness control
- 3 24 V power supply
- 4 A/B Port A/B for spot/ringlight
- 10 DSE connectors
- 12 segment switch
- 13 segment rotation
- (E) Backlight

# OPERATING INSTRUCTIONS FOR LED ILLUMINATION UNIT

Please read all the operating instructions before operating the LED illumination unit!

## Technical data

	LED ringlight 66/40	LED ringlight 66/80	LED ringlight 38/20
Diodes	40 white LEDs	80 white LEDs	20 white LEDs
Color temperature	5,000 K	5,000 K	5,000 K
Dimensions	outer: Ø 94 x 25 inner: Ø 66 mm	outer: Ø 112 x 23.5 mm inner: Ø 66 mm	outer: Ø 60 x 25 mm inner: Ø 38 mm
Attachment	Circumferential clamp with spring and screw		
	LED spot	LED backlight	
Diodes	19 white LEDs	20 white LEDs	
Color temperature	5,000 K	5,000 K	
Dimension	Ø 34 mm x 45 mm		
Attachment	using goose neck		
Carrier plate		PVC black, Ø 120 mm x 24.5 mm	
Object holder		compound glass, Ø 95 mm	
Illuminated surface		Ø 55 mm	
Control unit	Standard	Segment (40 LED)	Segment (80 LED)
Supply		24 V DC +/- 5 %, hollow plug 5.5 x 2.1 mm	
	250 mA	250 mA	500 mA
Dimension	95.5 x 64 x 29 mm	95.5 x 64 x 29 mm	102 x 106 x 39 mm
Connectors for ringlight/spot	2 x 9-pin miniature DIN socket	2 x 9-pin miniature DIN socket	D-Sub (25 pin)
On/Off switch	On/Off	with segment selector	with segment selector
Button for segment rotation	–	1	2
Potentiometer wheel		3–30 mA, with index holes	
ESD connection		4 mm banana connector	
Operating ambient temperature		10–40 °C (*)	
Rel. humidity		30–70 %	
Mark of conformity		CE	
Power supply			
Plug connector		3-pin IEC 320	
Voltage range		100–240 V	
Mains frequency		50–60 Hz	
Power rating max.		15 VA	
Protection class		I	
Mark of conformity		CE, UL, CSA	

(\*) LEDs are suitable for use at temperatures up to 40 °C. However, increasing temperature and humidity results in a reduction in lifetime. For optimal use, the ambient temperature should be 25 °C (or lower) and the relative humidity 50 % (or lower).

 CE Mark of conformity: Certifies that the control unit conforms to the EMC directive 89/336/EEC

## DESCRIPTION

This new type of LED illumination device (Light Emitting Diode) has been developed for stereomicroscopy illumination purposes in industrial and laboratory settings.

Please use the diagrams on the left-hand cover of this operating manual. In purchasing this product, you have acquired a high-quality LED illumination unit for all applications requiring intensive light with minimal heat generation. LEDs are employed as the illuminant.

The following features are characteristic of this type of illumination:

- ▼ Color temperature 5,000 K (daylight!)
- ▼ Ripple and flicker free
- ▼ Long LED lifetime
- ▼ Quiet operation, free of vibration
- ▼ Extremely compact design
- ▼ Battery operation possible

The modular design of the LED illumination device in combination with the

- ▼ Ringlight – provides homogeneous shadow-free illumination
- ▼ Spot – provides directional illumination
- ▼ Backlight – provides homogeneous illumination which passes the object

The illumination system consists of:

- ▼ (Desktop) power-supply
- ▼ Control unit and
- ▼ LED unit (1 ringlight 66/40 or 1 or 2 spots or backlight or ringlight 38/20)

The **power-supply** (C) voltage range extends from 100–240 V and features a 3-pin plug connector.

The compact **control unit** (D) has an On/Off switch (1), a continuously adjustable brightness control (2) with a „turbo range“ and two connector plugs (4) for connecting the LED illumination units.

A 24 V battery or an accumulator can be attached to the control unit via a special cable by means of the hollow plug connector (3).

The **LED ringlight** (40 LEDs) (A) allows surfaces to be illuminated exceptionally brightly and homogeneously and fits all microscopes with an objective diameter of 66 mm. Adapters are available for use with other microscopes. The ringlight is attached by means of a spring and screw, thereby protecting the objective.

The **LED spot** (19 LEDs) (B) can be attached directly to the microscope using a goose neck (6) or can be positioned as required using the goose neck and a base plate.

A wide range of accessories allows this LED illumination to be used for many purposes.

The **LED backlight** (20 LEDs) (E) can be used instead of the microscope's object holder to illuminate transparent objects from the bottom.

# SAFETY INFORMATION

Two danger symbols are used throughout this document.

Caution: Risk of electrical shock

Caution: Dangerous area.  
Warning: Refer to accompanying documentation.

This section contains safety information which must be observed strictly when using this device.

IT IS IN YOUR OWN INTERESTS TO PAY ATTENTION TO ALL WARNINGS on the unit and in this manual.

## Target audience

People working with the device must read the sections which are relevant to their work. This applies especially to the chapter entitled "Safety information".

## Duty of observation with respect to product

The operator must report all operational irregularities or changes to components which are relevant to device safety immediately to the responsible supervisor or the manufacturer.

## Location of operating instructions

We recommend storing these operating instructions near the device to ensure quick access by operating staff.

## Legal provisions

National and local safety and accident prevention regulations which are in force must be strictly observed in addition to the operational guidelines issued by the operating entity.

**Do not make any technical modifications to the device under any circumstances!**

Refer to specifications in "Technical Data" section for binding operational limits.

DO NOT USE this unit near water or in any area with excessive moisture.  
WARNING: In order to prevent electric shock, do not expose this appliance to rain or high humidity.

NEVER SPILL LIQUID ON THE UNIT OR INSERT OBJECTS INTO THE UNIT! This could result in electric shock or damage to the unit.

DO NOT PLACE FLAMMABLE MATERIALS on or near the unit at any time. Keep unit AWAY from other sources of HEAT. The device has not been approved for operation in areas subject to explosion hazards!

The device may only be operated using the mains voltage indicated.

NEVER OPEN THE APPLIANCE OR ANY COMPONENTS, unless instructed expressly to do so by these instructions.

Never look directly at the LEDs when switched-on, otherwise you risk eye injuries.

## Cleaning

Disconnect unit from the mains power supply before cleaning and only clean with a damp cloth. Never use combustible or flammable liquids. If fluids accidentally enter the device, unplug the mains cable and let the unit dry thoroughly before using again.

## Spare parts

Use only original spare parts. If this is not done, it can lead to personal injury and material damage. Refer to "Spare Parts" section for component codes.

## Responsibility

As the ultimate legal entity, the operating institution is responsible for ensuring the proper use of the device and for providing other operators with the necessary information and it specifies the competencies required to operate the device.

DISCONNECT THE MAINS CABLE when the unit is not being used for an extended period of time.

ONLY USE THE ORIGINAL MAINS CABLE. Route cable so that it cannot be jammed or severed.

FOR REASONS OF SAFETY only use the grounded 3-pin plug.

Repairs which are not described in this document must only be carried out by authorized workshops!

The manufacturer is not liable for any damage resulting from a failure to comply with the above instructions!

## Current technology

Photonic LED illumination devices constitute state-of-the-art technology and employ recognized safety standards.

This instrument was designed and built in accordance with the following regulations and standards:

- ▼ EEC/89/336: Electromagnetic compatibility (CE)
- ▼ ISO 9001/2: Quality systems for development and production

## INSTALLATION AND CONNECTION

Remove the individual system components from their packaging and place them onto a horizontal surface.

To obtain optimal functionality, the LED illumination system must be positioned correctly.

Observe the following criteria:

- ▼ No high atmospheric humidity (cf. Technical data)
- ▼ No high room temperature (cf. Technical data)

The device may only be operated at the rated mains voltage. Only connect the unit to grounded sockets.

Keep air vents free!

### **The control unit must only be operated with the power-supply provided!**

Connect one end of the power-supply to the control unit (3) (24 V DC input) by means of the hollow plug connector (9) and the other end to the power mains (100–240 V AC, 50–60 Hz) by means of the mains cable (8). If necessary, use the ESD connector (10) for potential equalisation.

Plug the LED lighting into the control unit. Only use Port A of the control unit for the LED ringlight 66/40 (and all future products with 40 LEDs). LED spot, backlight and ringlight 38/20 can be used on Port A and/or B.

The following possibilities are available with the two Ports A and B (4) of the control unit:

- ▼ 1 LED ringlight 66/40 on Port A
- ▼ 1 or 2 LED spots or LED backlights or LED ringlights 38/20 on Port A and B

The LED ringlight can now be mounted directly onto the objective of all microscopes with 66 mm diameter by means of a clamping ring which does not damage the surface of the objective. For objectives with any other diameter there is a wide range of adapters available.

For mounting the LED spot, "goosenecks" (6) of different lengths and threads are available, allowing the LED spot to be pointed in any desired direction. The one end (with the O-ring) of the gooseneck is inserted into the opening of the LED spot. The other end can be screwed directly onto the microscope or optionally into a base plate.

## OPERATION


The LED illumination unit can be switched on and off (Position I: On; Position 0: Off) using the switch located on the control unit (1).

The potentiometer wheel (2) is used to control the brightness, with the current being adjustable from 3 mA to 30 mA. Please note that the current level influences not only the brightness, but also the service life of the LEDs.

The approximate potentiometer wheel setting can be read from the scale on the wheel.

View of standard control unit

View of control unit with segment control

The following segments can also be selected with the on/off switch:

- ▼ Semicircle
- ▼ Quarter circle
- ▼ 2 opposing quarter circles

This segment selection can be rotated using (13).

Do not open the unit or its individual components.

### Lifetime

LED chips have a typical lifetime of about 100.000 hours. The actual service life depends to a great extent upon the ambient temperature and the operating current of the LEDs. A higher operating current produces not only more light but also more heat, which reduces the life span of the LEDs.

As a rule, the brightness of LEDs also decreases over the life cycle. After about 7000 hours at the highest possible luminance – i.e. at the potentiometer wheel "turbo" setting – the brightness will drop by 50 %. On the other hand, the same value will be reached in the 10 mA position only after 20,000 hours service. This figure increases substantially with a setting below 10 mA.


To increase the lifetime of the LEDs, the brightness should be adjusted to the actual level required at any point in time.

Faulty LEDs can be replaced. Contact your dealer or the manufacturer with respect to this matter.



## SPARE PARTS AND ACCESSORIES

Component description	Article number
Power supply 100–240 V AC, 24 V DC/15 W	ENSS-15-001
Control unit (max. 40 LED)	619-30-001
Control unit with segment control (max. 40 LED)	619-30-021
Control unit with segment control (max. 80 LED) for ringlight 66/80 only	619-30-022
LED ringlight 66/40 (40 LEDs)	619-20-002
LED ringlight 66/80 (80 LEDs)	619-20-035
LED ringlight 38/20 (20 LEDs)	619-20-006
LED spot (19 LEDs)	619-20-011
LED backlight (20 LEDs)	619-20-032
Gooseneck for spot M5 – 85 mm	619-30-003
Gooseneck for spot M5 – 200 mm	619-30-004
Gooseneck for spot M5 – 300 mm	619-30-005
Blue filter for spot	619-30-010
Green filter for spot	619-30-011
Yellow filter for spot	619-30-012
Red filter for spot	619-30-013
Base plate with M5 adapter for gooseneck	619-30-007
Polarization filter set for LED ringlight 66/40	619-20-021
Clamping ring for LED ringlight 66/40 + 66/80	595-32-013
Clamping ring screw for ringlight	619-30-008
LED ringlight adapter for:	
Objective 66/58 mm	595-32-258
Olympus SZ III, SZ 60, SZ 1145	598-30-002
Olympus SD, SF, SZ with lens attachment	598-30-017
Olympus SZH, SZX9, SZX12, Nikon SMZ 800, SMZ/U, 10A, 1000, 1500	598-32-043
Nikon (except for SMZ-1)	598-30-004
Olympus SD, SF, SZ without lens attachment, Nikon SMZ-1, 1B, 2B, 2T	598-32-049
Nikon SMZ 645 + SMZ 660	598-32-050
Universal adapter for diameters from 38 mm to 58 mm	598-30-022
Adapter DM38 / DM36 mm	619-32-056
Mains cable, 3-pin IEC320 for:	
EU	ELN-3.001
CH	ELN-3.002
IT	ELN-3.009
UK	ELN-3.004
USA/CAN/J	ELN-3.003
AUS	ELN-3.006






The design of the product is under constant review, and although every effort is made to keep this manual up to date, we reserve the right to change specifications at any time without prior notice. Not responsible for typographical and misprint errors.

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