

MICROSCOPE COMPONENTS GUIDE



Choosing The Ideal UIS2 Optics Components For Your Equipment





The wide range of Olympus components introduced here allows users in such diverse fields as research, inspection and production to take advantage of the quality, flexibility and outstanding optical performance of the UIS2 Optical System.

That's why installing Olympus microscope components is, quite simply, the right choice for your equipment.

CONTENTS

	OME TO UIS2/UIS OPTICS			
SYSTI	EM DIAGRAM		5 -	<u> </u>
UIS2/	UIS OBJECTIVE LENSES			
	M Plan SemiApochromat	MPLFLN series		8
	Long WD M Plan SemiApochroma	t LMPLFLN series		9
	M Plan Achromat	MPLN series		- 10
	LCD Long WD M Plan SemiApochr	omat LCPLFLN-LCD series		11
	M Plan Apochromat	MPlanApo series		12
	Super Long WD M Plan Achromat	SLMPlan series		12
	IR Long WD M Plan SemiApochror	mat		
		LMPlan-IR series		
	IR M Plan SemiApochromat M Plan SemiApochromat BD	MPlan-IR MPLFLN-BD series		
	M Plan SemiApochromat BDP	MPLFLN-BDP series		1: 1:
	Long WD M Plan SemiApochroma			
		LMPLFLN-BD series		
	M Plan Achromat BD	MPLN-BD series		17
	M Plan Apochromat BD	MPlanApo-BD		18
MICR	OSCOPE SYSTEM BXFM		19 —	- 22
	BXFM frame	BXFM-F		
	BXFM	BXFM-F+BXFM-ILH+BXFM-ILHSPU		
	BXFM-S	BXFM-F+BXFM-ILHS		
	Universal stand type 2	SZ2-STU2		
	Compact stand Large stand	SZ-STL		
	3			
	AINIATION LINITS		າາ	21
ILLUIV			23 —	- 2:
	Reflected light illuminator for BF/E	лг BX-RLA2		- 23
	Universal reflected light illuminate	rBX-URA2		- 23
	Reflected light illuminators for BF	BX-KMA/BX-KMA-ESD		- 24
	Reflected light illuminator for BF			- 25
LAMP				
	75W xenon apo lamp housing	U-LH75XEAPO		
	100W mercury apo lamp housing	U-LH100HGAPO U-LH100HG		
	100W mercury lamp housing 100W halogen lamp housings	U-LH100-3/U-LH100IR/U-LH100L-3		
	External power supply	TH4-100/200		
	Hand switch	TH4-HS		
	Extension cord	U-RMT		
	DF converter for BX-URA2	U-RCV		- 27
	Fiber adapter for reflected light ob			
	Transmitted light guide adapter	U-LGAD SZX-TLGAD		
	Transmitted light guide adapter Light source	LG-PS2		
	Light guide	LG-SF		
	Double lamp house adapter	U-DULHA		
OBSE	RVATION TUBES		29 —	- 30
JUJE	Widefield hinocular tube	II_TD30_2		_ 20
JJJL	Widefield binocular tube Widefield binocular tube for IR	U-TR30-2		
JUJE	Widefield binocular tube for IR	U-TR30IR		
JJJL		U-TR30IR		- 29
JJ02	Widefield binocular tube for IR	U-TR30IR ube U-ETR-4 U-TLU		- 29 - 29 - 29
JJ02	Widefield binocular tube for IR Widefield erect image trinocular to Single port tube with lens Single port tube with lens for IR	U-TR30IR ube U-ETR-4 U-TLU U-TLUIR		- 29 - 29 - 29
5552	Widefield binocular tube for IR Widefield erect image trinocular tu Single port tube with lens Single port tube with lens for IR Super widefield trinocular tube	U-TR30IR ube U-ETR-4 U-TLU U-TLUIR U-SWTR-3		- 29 - 29 - 29
5552	Widefield binocular tube for IR Widefield erect image trinocular to Single port tube with lens Single port tube with lens for IR	U-TR30IR		- 29 - 29 - 29 - 30
JUSE	Widefield binocular tube for IR Widefield erect image trinocular tu Single port tube with lens Single port tube with lens for IR Super widefield trinocular tube Super widefield erect image trinocular	U-TR30IR		- 29 - 29 - 29 - 30
SSSE	Widefield binocular tube for IR Widefield erect image trinocular tu Single port tube with lens Single port tube with lens for IR Super widefield trinocular tube	U-TR30IR		- 29 - 29 - 29 - 30
	Widefield binocular tube for IR Widefield erect image trinocular tu Single port tube with lens Single port tube with lens for IR Super widefield trinocular tube Super widefield erect image trinocular	U-TR30IR ube U-ETR-4 U-TLU U-SWTR-3 U-SWTR-3 U-SWETR U-SWETR		- 29 - 29 - 29 - 30
	Widefield binocular tube for IR Widefield erect image trinocular tu Single port tube with lens Single port tube with lens for IR Super widefield trinocular tube Super widefield erect image trinoc Super widefield erect image tilting	U-TR30IR ube U-ETR-4 U-TLU U-SWTR-3 U-SWTR-3 U-SWETR U-SWETR		- 29 - 29 - 29 - 30 - 30
	Widefield binocular tube for IR Widefield erect image trinocular tu Single port tube with lens Single port tube with lens for IR Super widefield trinocular tube Super widefield erect image trinoc Super widefield erect image tilting	U-TR30IR	31 —	- 29 - 29 - 29 - 30 - 30 - 31 - 31
	Widefield binocular tube for IR Widefield erect image trinocular tu Single port tube with lens Single port tube with lens for IR Super widefield trinocular tube Super widefield erect image trinoc Super widefield erect image tilting RMEDIATE TUBES & ACCESS Magnification changer Magnification changer 2x	U-TR30IR	31 —	- 29 - 29 - 29 - 30 - 30 - 31 - 31
	Widefield binocular tube for IR Widefield erect image trinocular tu Single port tube with lens Single port tube with lens for IR Super widefield trinocular tube Super widefield erect image trinoc Super widefield erect image tillting RMEDIATE TUBES & ACCESS Magnification changer	U-TR30IR	31 —	- 29 - 29 - 29 - 30 - 30 - 30 - 31 - 31
	Widefield binocular tube for IR Widefield erect image trinocular tu Single port tube with lens Single port tube with lens for IR Super widefield trinocular tube Super widefield erect image trinoc Super widefield erect image tilting RMEDIATE TUBES & ACCESS Magnification changer 2x Trinocular intermediate attachmen	U-TR30IR	31 —	- 29 - 29 - 29 - 30 - 30 - 30 - 31 - 31 - 31
	Widefield binocular tube for IR Widefield erect image trinocular tu Single port tube with lens Single port tube with lens for IR Super widefield trinocular tube Super widefield erect image trinoc Super widefield erect image tilling RMEDIATE TUBES & ACCESS Magnification changer Magnification changer 2x Trinocular intermediate attachmen Dual port	U-TR30IR	31 —	- 29 - 29 - 29 - 30 - 30 - 31 - 31 - 31 - 31
	Widefield binocular tube for IR Widefield erect image trinocular tu Single port tube with lens Single port tube with lens for IR Super widefield trinocular tube Super widefield erect image trinoc Super widefield erect image tilting RMEDIATE TUBES & ACCESS Magnification changer 2x Trinocular intermediate attachmen	U-TR30IR	31 —	- 29 - 29 - 29 - 30 - 30 - 30 - 31 - 31 - 32 - 32 - 32

EYEPIECES/FILAR MICROMETER	EYEPIECE	3
Widefield eyepieces	WHN10x, WHN10x-H, CROSSWHN10x, WH15x	
Super widefield eyepieces	SWH10x-H, MICROSWH10x, CROSSSWH10x	
Filar micrometer eyepiece	U-OSM	
REVOLVING NOSEPIECES	34 -	
Quintuple revolving nosepiece Sextuple revolving nosepiece with	U-5RE-2 n slider slot for DIC U-D6RE	
	o slider slot for DIC with ESD treatment U-D6RE-ESD	
Septuple revolving nosepiece with		
Centerable quadruple revolving no		
Centerable sextuple revolving nos	epiece with slider slot for DIC U-P6RE	3
Quintuple revolving nosepiece for	BF/DF U-5BDRE	3
Quintuple revolving nosepiece for	U-D5BDRE	3
Sextuple revolving nosepiece for I	U-D6BDRE	3
Centerable quintuple revolving no	sepiece U-P5BDRE	3
Adapter to mount BF objectives	BD-M-AD	3
VIDEO CAMERA ADAPTERS	36 -	— — 3
C-mount video camera ports	U-TV0.25xC, U-TV0.35xC-2, U-TV0.5xC-3, U-TV0.63xC	3
Video camera mount adapters		
Video camera port	U-TV1x-2	3
MOTORIZED LINUTO	20	_
	38 - uminator+motorized Nomarski DIC sextuple	
revolving nosepiece+100W halog	en lamp housing	
Motorized universal reflected ligh	BX-RLAA+U-D6REMC+U-LH100-3 t illuminator BX-RFAA	
Motorized quintuple BD revolving		
Motorized sextuple revolving nose	epiece with slider slot for DIC U-D6REMC	?
Motorized centerable quintuple re	volving nosepiece with slider slot for DIC U-P5REMC	
Control unit	BX-UCB	3
Hand switch Control box for motorized nosepie	U-HSTR2	3
control box for motorized nosepie	BX-REMCB	3
AC adapter for BX-REMCB	U-ACAD4515	3
Active auto focus unit	U-AFA1M	
Motorized reflected filter wheel	U-FWR	4
Motorized illumination with power	BXFMA-F	4
DEEP UI TRAVIOI ET OBSERVATIO	ON SYSTEM	
UV248 compatible intermediate to		
2.2.12 11panbio intorniodiato te	U-UVF248IM	4
UV quartz light guide	U-UVF2FV/5FB	4
UV248 compatible light source bo	x + Mercury Xenon lamp housing U-UVF248LB+U-LH80HBXE	4
OPTICAL TERMINOLOGY	43	
	40	


WELCOME TO UIS2/UIS OPTICS

UIS2/UIS:

The System That Maximizes The Advantage Of Infinity-Corrected Optics

What's infinity-corrected optics?

UIS2/UIS optics is an infinity-corrected optical system – in other words, a system in which light passes from the specimen through the objective lens without forming an image along the way. Instead, it travels in the form of infinity parallel rays to the tube lens. The tube lens is where the intermediate image is formed, whereas in finite-corrected optics, this is done by the objective lens.



Advantages of infinity-corrected optics


This system, known as "infinity-corrected optics", offers a number of advantages:

- There is no change in magnification even when the distance between the objective lens and tube lens is altered.
- With the total magnification remaining constant, there is no image aberration – even when prisms or sliders are interposed between the objective lens and the tube lens.

As thousands of users have found by experience, these advantages are crucial to composing the ideal microscope optical system. What's more, it is even possible to freely insert or remove intermediate attachments in the parallel rays of

light between the objective lens and tube lens, allowing the creation of user-specific or task-specific optical systems. To establish real flexibility with such a system, it is necessary to eliminate the occurrence of coma aberration.


*In UIS2/UIS objective lenses, the parfocal distance is designed at 45mm and the focal length of the tube lens is 180mm.



Basic dimensions of UIS2/UIS optical system

The UIS2/UIS optical system optimally corrects aberration with a dedicated telan lens and an eyepiece so that the coma aberration and flatness are not degraded even when the telan lens exit pupil position is changed by changing the objective lens and telan distance. This makes it possible to use a distance of 50mm to 170mm from objective lens mounting position to the single port tube with lens.

*Coma aberration: refer to the optical terminology at the end of this document.



WELCOME TO UIS2/UIS OPTICS

Features of UIS2 objective lenses

UIS2 objective lenses ensure compatibility (screw diameter, optical performance) with the UIS optical system and have the following features compared to conventional objective lenses.

1. Wavefront aberration control

The Olympus UIS2 objective lenses set a new standard, with wavefront aberration control in addition to common performance standards of N.A. and W.D. Olympus challenges farther highest order optics which has not been fulfilled by the conventional standards. We offer excellent performance objective lenses by minimizing the aberrations that lower resolution. *Wave front aberration: refer to the optical terminology at the end of this document.

2. Objective lenses with excellent image parcentricity

High power SemiApochromatic UIS2 objective lenses make the centration tolerance between objective lenses on the microscope nosepiece keep the image within the enter of the field of view even with digital cameras. (50x or higher power in both MPLFLN and LMPLFLN series)

3. Improvement of color reproducibility

UIS2 objective lenses realize natural color reproduction without any chromatic shifts using stringently selected high transmittance glass and advanced coating technology that provides high transmittance which is flat over an ultrawide band wavelength. In addition, since the total optical system, including the tube lens is designed to reproduce a natural color, clear images faithful to the specimen are obtained even with digital imaging.

4. Lightening

Weight has been reduced to approximately 2/3 that of conventional products by using an aluminum objective lens barrel cover. This has the effect of lightening the load on the devices at objective lens up/down, suppressing vibrations by lowering the inertial force at objective lens switching, etc. (MPLFLN series, LMPLFLN series)


5. Adoption of eco-lens

The glass materials of UIS2 objective lenses are all lead- and cadmium-free eco-glass.


Based on our conviction that the UIS2/UIS system is the best way to maximize the advantages of infinity-corrected optical systems, we confidently recommend the UIS2/UIS-featured Olympus microscope units for all your high-precision needs in research, inspection and production equipment.


* Refer to the Olympus home page for detailed objective lenses specifications.

SYSTEM DIAGRAM




SYSTEM DIAGRAM





■ Meaning of abbreviations shown on objective lens



■ Objective lens series list

	Series	Magnification	BF	DF	DIC ⁻¹	POL	FL	F.N. (Field Number)	Remarks
UIS2	MPLFLN	1.25/2.5	0					1.25x: 22 / 2.5x: 26.5	Use together with polarizer and analyzer
									recommended
		5/10/20/50/100	0		OU	0	0,5	26.5	
	LMPLFLN	5/10/20/50/100	0		OL	0	0	26.5	
	MPLN	5/10/20/50/100	0					22	
	LCPLFLN-LCD	20/50/100	0		OL			26.5	For LCD
UIS	MPlanApo	20/50/100	0		O.3 N	0		26.5	
	SLMPlan	20/50	0					26.5	
	LMPlanIR/MPlanIR	5/10/20/50/100 *4	0					22	For near-IR observation
UIS2	MPLFLN-BD	5/10/20/50/100/150	0	0	OU	0	O*2	26.5	
	MPLFLN-BDP	5/10/20/50/100	0	0	●U	•	O*2	26.5	
	LMPLFLN-BD	5/10/20/50/100	0	0	OL	0	0	26.5	
	MPLN-BD	5/10/20/50/100	0	0				22	
UIS	MPlanApo-BD	100	0	0	OU	0		26.5	

^{**} DIC prism U-DICR: UM/LM position, U-DICRHC: LM position fixed, U-DICRH: UM position fixed. *2 5~20x: U excitation also possible

DIC: Differential Interference Contrast POL: Polarized light FL: Fluorescence

Cover glass thickness (no cover)

■ Features of objective lens series

MPLFLN series: M Plan SemiApochromat — P 8

Plan SemiApochromat objective lenses, giving high-level correction for chromatic aberration. The lineup consists of 7 objective lenses ranging from 1.25x to 100x, and secures a W.D. of 1mm or longer. Since the exit pupil position of the 5x-100x objective lenses is standardized, the position of the DIC prism does not have to be switched when changing the magnification. For ultra low magnifications (1.25x, 2.5x), use together with analyzer and polarizer of the reflected light illuminator.

• LMPLFLN series: Long WD M Plan SemiApochromat — P 9

Long working distance Plan SemiApochromat objective lenses, giving high-level correction for chromatic aberration. Suitable with samples having a height difference and in preventing collision, as the working distance is long. Also, since the exit pupil position of the 5x-100x objective lenses is standardized, the position of the DIC prism does not have to be switched when changing the magnification.

MPLN series: M Plan Achromat — P 10

Plan Achromat objective lenses providing excellent image flatness up to F.N. 22.

• LCPLFLN-LCD series: LCD Long WD M Plan SemiApochromat — P 11

Perfect objective lens series for observation of LCD panels and other samples through a glass substrate. Aberration correction matched to the glass thickness is accomplished using a correction ring.

• MPlanApo series: M Plan Apochromat — P 12

Highest class Plan Apochromat objective lenses that maximize performance in brightfield observation. All aberrations are corrected at the highest level, while providing high N.A.

• SLMPlan series: Super Long WD M Plan Achromat — P 12

Plan Achromat objective lenses with high magnification and super long working distance. Two magnifications, 20x and 50x are available. For 5x or 10x objective lenses, select from the LMPLFLN Series.

LMPlan-IR series: IR Long WD M Plan SemiApochromat — P 13 MPlan-IR: IR M Plan SemiApochromat — P 13

IR objective lenses which compensate for aberrations from visible to near infrared light. Ideal for the observations of semiconductor interiors and the back surface of a chip package as well as CSP bump inspection.

MPLFLN-BD series: M Plan SemiApochromat BD — P 14

Field Number

Plan SemiApochromat objective lenses, giving high-level correction for chromatic aberration. The series secures a W.D. of 1mm or longer. Since the exit pupil position of the 5x-150x objective lenses is standardized, the position of the DIC prism does not have to be switched when changing the magnification.

MPLFLN-BDP series: M Plan SemiApochromat BDP — P 15

Plan SemiApochromat objective lenses, giving high-level correction for chromatic aberration. The series secures a W.D. of 1mm or longer. Since the exit pupil position of the 5x-100x objective lenses is standardized, the position of the DIC prism does not have to be switched when changing the magnification. The BDP series optimizing brightfield/darkfield and polarized light characteristics is perfect for Nomarski DIC and polarized light observations.

■ LMPLFLN-BD series: Long WD M Plan SemiApochromat BD — P 16

Long working distance Plan SemiApochromat objective lenses, giving high-level correction for chromatic aberration. Suitable with samples having a height difference and in preventing collision, as the working distance is long. Also, since the exit pupil position of the 5x-100x objective lenses is standardized, the position of the DIC prism does not have to be switched when changing the magnification.

MPLN-BD series: M Plan Achromat BD — P 17

Plan Achromat objective lenses providing excellent image flatness up to F.N. 22.

MPlanApo BD: Plan Apochromat BD — P 18


Highest class Plan Apochromat objective lens that maximize performance in brightfield and darkfield observations. All aberrations are corrected at the highest level, while providing bigh N A

^{*3 50}x: DIC observation not applicable
*4 MPlanIR: available 100x only ○: Responds ●: Optimally responds BF: Brightfield DF: Darkfield

M Plan SemiApochromat

MPLFLN series

Plan SemiApochromat objective lenses, giving high-level correction for chromatic aberration. The lineup consists of 7 objective lenses ranging from 1.25x to 100x, and secures a W.D. of 1mm or longer. Since the exit pupil position of the 5x-100x objective lenses is standardized, the position of the DIC prism does not have to be switched when changing the magnification. For ultra low magnifications (1.25x, 2.5x), use together with analyzer and polarizer of the reflected light illuminator.




	UIS2	objective lenses			Widefield eyepiece WHN10x Field Number 22			Super widefield eyepiece SWH10x Field Number 26.5		
Objective lens (magnification)	Numerical Aperture	Working distance (mm)	Focal distance f (mm)	Weight (g)	Total magnifications	Practical field of view (mm)	Depth of focus (µm)	Total magnifications	Practical field of view (mm)	Depth of focus (µm)
MPLFLN 1.25x*	0.04	3.5	145	122	12.5	17.6	870	_	_	_
MPLFLN 2.5x*	0.08	10.7	72	106	25	8.8	220	25	10.6	220
MPLFLN 5x	0.15	20.0	36	51.5	50	4.4	59	50	5.3	59
MPLFLN 10x	0.30	11.0	18	68.1	100	2.2	15	100	2.7	15
MPLFLN 20x	0.45	3.1	9	70.4	200	1.1	5.2	200	1.3	5.1
MPLFLN 50x	0.80	1.0	3.6	89.9	500	0.44	1.3	500	0.53	1.3
MPLFLN 100x	0.90	1.0	1.8	90.9	1000	0.22	0.73	1000	0.27	0.73


Long WD M Plan SemiApochromat


(WD: Working Distance)

LMPLFLN series


Long working distance Plan SemiApochromat objective lenses, giving high-level correction for chromatic aberration. Suitable with samples having a height difference and in preventing collision, as the working distance is long. Also, since the exit pupil position of the 5x-100x objective lenses is standardized, the position of the DIC prism does not have to be switched when changing the magnification.







LMPLFLN50x




Unit: mm


	UIS2	objective lenses		Widefield eyepiece WHN10x Field Number 22			Super widefield eyepiece SWH10x Field Number 26.5			
Objective lens (magnification)	Numerical Aperture	Working distance (mm)	Focal distance f (mm)		Total magnifications	Practical field of view (mm)	Depth of focus (µm)	Total magnifications	Practical field of view (mm)	Depth of focus (µm)
LMPLFLN 5x	0.13	22.5	36	50	50	4.4	70	50	5.3	70
LMPLFLN 10x	0.25	21.0	18	54	100	2.2	18	100	2.7	18
LMPLFLN 20x	0.40	12.0	9	73	200	1.1	6.1	200	1.3	6.1
LMPLFLN 50x	0.50	10.6	3.6	77	500	0.44	2.5	500	0.53	2.5
LMPLFLN 100x	MPLFLN 100x 0.80 3.4 1.8 94						0.87	1000	0.27	0.87

M Plan Achromat

MPLN series

Plan Achromat objective lenses providing excellent image flatness up to F.N. 22.






Unit: mm

				Widefield eyepiece WHN10x Field Number 22				
	Objective lens Numerical (magnification) Aperture		Working distance (mm)	Focal distance f (mm)	Weight (g)	Total magnifications	Practical field of view (mm)	Depth of focus (µm)
MPLN	5x	0.10	20.0	36	64	50	4.4	98
MPLN	10x	0.25	10.6	18	80	100	2.2	18
MPLN	20x	0.40	1.3	9	111	200	1.1	6.1
MPLN	50x	0.75	0.38	3.6	113	500	0.44	1.4
MPLN	100x	0.90	0.21	1.8	116	1000	0.22	0.73

LCD Long WD M Plan SemiApochromat

LCPLFLN-LCD series

Perfect objective lens series for observation of LCD panels and other samples through a glass substrate. Aberration correction matched to the glass thickness is accomplished using a correction ring.



^{*} Value at glass thickness 0.7mm observation

Objective lens	LCPLFLN20xLCD			L	CPLFLN50xL0	CD	LCPLFLN100xLCD			
Corresponding glass thickness (mm)	0-1.2				0-1.2		0-0.7			
Correction ring indication	0	0.7	1.2	0	0.7	1.2	0	0.5	0.7	
Working distance (mm)	8.3	7.8	7.4	3.0	2.5	2.2	1.2	0.98	0.9	
Correction system	Correction ring		Correction ring			Correction ring				


	UIS2	objective lenses		Widefield eyepiece WHN10x Field Number 22			Super widefield eyepiece SWH10x Field Number 26.5			
Objective lens (magnification)	Numerical Aperture	Total magnifications	Practical field of view (mm)		Total magnifications	Practical field of view (mm)	Depth of focus (µm)			
LCPLFLN 20xLCD***	0.45	7.8	9	146	200	1.1	5.2	200	1.3	5.2
LCPLFLN 50xLCD***	0.70	2.5	3.6	170	500	0.44	1.6	500	0.53	1.6
LCPLFLN 100xLCD***	0.85	0.9	1.8	185	1000	0.22	0.79	1000	0.27	0.79


^{**}The figure shown here is the value when the correction ring indication is 0.7. *** To be available in the beginning of 2007


M Plan Apochromat

MPlanApo series

Highest class Plan Apochromat objective lens that maximize performance in brightfield observation. All aberrations are corrected at the highest level, while providing high N.A.







Unit: mm

	UIS o	bjective lenses	Widefield eyepiece WHN10x Super widefield eyepiece Field Number 22 Field Number 26.5							
Objective lens Numerical Working distance Focal distance Weight (magnification) Aperture (mm) f (mm) (g)					Total magnifications	Practical field of view (mm)		Total magnifications	Practical field of view (mm)	Depth of focus (µm)
MPlanApo 20x	0.60	0.9	9	150	200	1.1	3.7	200	1.3	3.7
MPlanApo 50x	0.95	0.3	3.6	150	500	0.44	1.0	500	0.53	1.0
MPlanApo 100x	0.95	0.35	1.8	150	1000	0.22	0.67	1000	0.27	0.67


Screw: W20.32x0.706 (0.8"x1/36")


Super Long WD M Plan Achromat

SLMPlan series

Plan Achromat objective lenses with high magnification and super long working distance.

Two magnifications, 20x and 50x are available. For 5x or 10x objective lenses, select from the LMPLFLN series.





Unit: mm


	UIS objective lenses						Widefield eyepiece WHN10x Super widefield eyepie Field Number 22 Field Number			
Objective lens Numerical Working distance Focal distance Weight (magnification) Aperture (mm) f (mm) (g)					Total magnifications	Practical field of view (mm)	Depth of focus (µm)	Total magnifications	Practical field of view (mm)	
SLMPlan 20x	0.35	21.0	9	73	200	1.1	7.2	200	1.3	7.2
SLMPlan 50x	0.45	15.0	3.6	91	500	0.44	2.9	500	0.53	2.9


IR Long WD M Plan SemiApochromat/IR M Plan SemiApochromat


LMPlan-IR series/MPlan-IR


IR objective lenses which compensate for aberrations from visible to near infrared light.


Ideal for the observations of semiconductor interiors and the back surface of a chip package as well as CSP bump inspection.














Unit: mm

			UIS objective lenses		Widefield eyepiece WHN10x Field Number 22				
	Objective lens Num (magnification) Ape		Working distance (mm)	Focal distance f (mm)	Weight (g)	Total magnifications	Practical field of view (mm)	Depth of focus (µm)	
LMPlan	5xIR	0.10	20.0	36	73	50	4.4	98	
LMPlan	10xIR	0.25	18.5	18	73	100	2.2	18	
LMPlan	20xIR	0.40	8.1	9	110	200	1.1	6.1	
LMPlan	50xIR	0.55	6.0	3.6	115	500	0.44	2.2	
LMPlan	100xIR	0.80	3.4	1.8	122	1000	0.22	0.87	
MPlan	100xIR	0.95	0.3	1.8	130	1000	0.22	0.67	


M Plan SemiApochromat BD

(BD:Brightfield/Darkfield)

MPLFLN-BD series

Plan SemiApochromat objective lenses, giving high-level correction for chromatic aberration.

The series secures a W.D. of 1mm or longer. Since the exit pupil position of the 5x-150x objective lenses is standardized, the position of the DIC prism does not have to be switched when changing the magnification.



Unit: mm


	UIS2	objective lenses			Widefield eyepiece WHN10x Super widefield eyep Field Number 22 Field Numbe					
Objective lens (magnification)	Numerical Aperture	Working distance (mm)	Focal distance f (mm)	Weight (g)	Total magnifications	Practical field of view (mm)	Depth of focus (µm)	Total magnifications	Practical field of view (mm)	Depth of focus (µm)
MPLFLN 5xBD	0.15	12.0	36	95.5	50	4.4	59	50	5.3	59
MPLFLN 10xBD	0.30	6.5	18	82.8	100	2.2	15	100	2.7	15
MPLFLN 20xBD	0.45	3.0	9	87.7	200	1.1	5.2	200	1.3	5.2
MPLFLN 50xBD	0.80	1.0	3.6	99.8	500	0.44	1.3	500	0.53	1.3
MPLFLN 100xBD	0.90	1.0	1.8	98.9	1000	0.22	0.73	1000	0.27	0.73
MPLFLN 150xBD	0.90	1.0	1.2	104.8	1500	0.15	0.6	1500	0.18	0.6

M Plan SemiApochromat BDP

(BDP:Brightfield/Darkfield/Polarizing)


MPLFLN-BDP series

Plan SemiApochromat objective lenses, giving high-level correction for chromatic aberration. The series secures a W.D. of 1mm or longer. Since the exit pupil position of the 5x-100x objective lenses is standardized, the position of the DIC prism does not have to be switched when changing the magnification. The BDP series optimizing brightfield/darkfield and polarized light characteristics is perfect for Nomarski DIC and polarized light observations



MPLFLN50xBDP

MPLFLN100xBDP




Unit: mm

	UIS2 objective lenses				Widefield eyepiece WHN10x Field Number 22			Super widefield eyepiece SWH10x Field Number 26.5		
Objective lens (magnification)	Numerical Aperture	Working distance (mm)	Focal distance f (mm)	Weight (g)	Total magnifications	Practical field of view (mm)	Depth of focus (µm)	Total magnifications	Practical field of view (mm)	Depth of focus (µm)
MPLFLN 5xBDP	0.15	12.0	36	95.5	50	4.4	59	50	5.3	59
MPLFLN 10xBDP	0.25	6.5	18	83.3	100	2.2	18	100	2.7	18
MPLFLN 20xBDP	0.40	3.0	9	88.5	200	1.1	6.1	200	1.3	6.1
MPLFLN 50xBDP	0.75	1.0	3.6	100.5	500	0.44	1.4	500	0.53	1.4
MPLFLN 100xBDP	0.90	1.0	1.8	101.5	1000	0.22	0.73	1000	0.27	0.73

Long WD M Plan SemiApochromat BD

LMPLFLN-BD series

Long working distance Plan SemiApochromat objective lenses, giving high-level correction for chromatic aberration. Suitable with samples having a height difference and in preventing collision, as the working distance is long. Also, since the exit pupil position of the 5x-100x objective lenses is standardized, the position of the DIC prism does not have to be switched when changing the magnification.




Unit: mm


	UIS2 objective lenses				Widefield eyepiece WHN10x Su Field Number 22			Super widefield eyepiece SWH10x Field Number 26.5		
Objective lens (magnification)	Numerical Aperture	Working distance (mm)	Focal distance f (mm)	Weight (g)	Total magnifications	Practical field of view (mm)	Depth of focus (µm)	Total magnifications	Practical field of view (mm)	Depth of focus (µm)
LMPLFLN 5xBD	0.13	15.0	36	81	50	4.4	70	50	5.3	70
LMPLFLN 10xBD	0.25	10.0	18	84	100	2.2	18	100	2.7	18
LMPLFLN 20xBD	0.40	12.0	9	86	200	1.1	6.1	200	1.3	6.1
LMPLFLN 50xBD	0.50	10.6	3.6	85	500	0.44	2.5	500	0.53	2.5
LMPLFLN 100xBD	0.80	3.3	1.8	102	1000	0.22	0.87	1000	0.27	0.87


M Plan Achromat BD

MPLN-BD series

Plan Achromat objective lenses providing excellent image flatness up to F.N. 22.









MPLN50xBD

MPLN100xBD





Unit: mm


		Wide	field eyepiece WH Field Number 22	N10x			
Objective lens (magnification)	Numerical Aperture	Working distance (mm)	Focal distance f (mm)	Weight (g)	Total magnifications	Practical field of view (mm)	Depth of focus (µm)
MPLN 5xBD	0.10	12.0	36	137	50	4.4	98
MPLN 10xBD	0.25	6.5	18	155	100	2.2	18
MPLN 20xBD	0.40	1.3	9	162	200	1.1	6.1
MPLN 50xBD	0.75	0.38	3.6	157	500	0.44	1.4
MPLN 100xBD	0.90	0.21	1.8	160	1000	0.22	0.73

M Plan Apochromat BD

MPlanApo-BD

Highest class Plan Apochromat objective lens that maximize performance in brightfield and darkfield observations. All aberrations are corrected at the highest level, while providing high N.A.

MPlanApo100xBD




Unit: mm

	UIS objective lenses					d eyepiece V eld Number 2			efield eyepied Id Number 2	
Objective lens (magnification)	3			Total magnifications	Practical field of view (mm)	Depth of focus (µm)	Total magnifications	Practical field of view (mm)	Depth of focus (µm)	
MPlanApo100xBD	0.9	0.31	1.8	180	1000	0.22	0.59	1000	0.27	0.59

BXFM frame

BXFM-F


Widely used system that allows use in combination with fiber illumination, motorized revolving nosepiece and telan lens unit. Can easily be integrated into other equipment. Attach to the equipment by rear bolt mounting screw or pillar mounting hole.




BXFM

BXFM-F+BXFM-ILH+BXFM-ILHSPU

Accommodates the reflected light brightfield/darkfield and fluorescence illuminators.




Weight: 3.2kg




BXFM-S

BXFM-F+BXFM-ILHS

Compact focusing unit suitable for building into existing equipment.






Stands

A wide variety of stands are available to suit different applications and purposes.


SZ2-STU2 Universal stand type 2




Major specifications

101	wajor specifications						
	Item	Specifications					
1	Diameter of focusing arm or fixing section of tube	ø32mm					
2	Vertical pole diameter	ø40mm					
3	Horizontal poles diameters	ø25mm					
		(both upper and lower poles)					
4	Stroke	Horizontal: 234mm,					
		Vertical: 205mm					
5	Movement range	Horizontal: 541 (435+106) mm max.					
		(Vertical pole —					
		BXFM-S optical axis)					
6	Maximum specimen weight	Forward: 10kg					
		(within 90-degree area)					
		Transverse direction: 6kg					
		Backward direction: 7kg					
		(at maximum stroke)					
7	Weight	30kg					

U-ST Compact stand



SZ-STL Large stand



Weight: 1.8kg Weight: 5kg " Unit: mm

 $^{^{\}star}$ The rotation angle of the horizontal arm can restrict to 90 degrees with stopper.

ILLUMINATION UNITS


Reflected light illuminator for BF/DF

BX-RLA2

ND filters are linked when exchanging between brightfield and darkfield.

Accessories

Unit name	Description	Weight (g)
U-25LBD	LBD filter slider	20
U-25IF550	IF550 filter slider	20
U-25ND6	ND filter	20
U-25ND25	ND filter	20
U-25FR	Frost filter slider	20
U-25L42	UV-cut filter	20
U-PO3	Polarizer slider for reflected light	71
U-POTP3	Polarizer slider for reflected light	71
	with tint plate	
U-AN360-3	360° rotatable analyzer slider	79
U-AN	Analyzer slider for reflected light	50
U-DICR	DIC slider for reflected light	130
U-DICRH	DIC slider for reflected light	130
	(high resolution type)	
U-DICRHC	DIC slider for reflected light	130
	(high contrast type)	



Weight: 3.4kg

Universal reflected light illuminator

BX-URA2

Suitable for observations ranging from brightfield to fluorescence.

Six mirror units can be attached to this reflected light illuminator simultaneously.

Accessories

Unit name	Description	Weight (g)
U-25LBD	LBD filter slider	20
U-25IF550	IF550 filter slider	20
U-25ND6	ND filter	20
U-25ND25	ND filter	20
U-25FR	Frost filter slider	20
U-25L42	UV-cut filter	20
U-PO3	Polarizer slider for reflected light	71
U-POTP3	Polarizer slider for reflected light	71
	with tint plate	
U-AN360-3	360° rotatable analyzer slider	79
U-AN	Analyzer slider for reflected light	50
U-DICR	DIC slider for reflected light	130
U-DICRH	DIC slider for reflected light	130
	(high resolution type)	
U-DICRHC	DIC slider for reflected light	130
	(high contrast type)	
U-MBF3	Mirror unit for reflected brightfield	80
U-MDF3*	Mirror unit for reflected darkfield	80
U-MDIC3	Mirror unit for reflected DIC	80
U-MBFL3	Mirror unit for reflected brightfield,	80
	for high intensity light source	
U-MWUS3	Fluorescence mirror unit for	80
	reflected (U excitation)	
U-MWBS3	Fluorescence mirror unit for	80
	reflected (B excitation)	
U-MWGS3	Fluorescence mirror unit for	80
	reflected (G excitation)	

 $[\]oplus$ **#** \oplus Ф (152)367

Weight: 3.8kg

^{*} U-RCV (DF converter for BX-URA2) is needed with darkfield observation.

ILLUMINATION UNITS


Reflected light illuminators for BF


BX-KMA/BX-KMA-ESD

Enables brightfield, Nomarski DIC and simple polarizing observations. ESD model is also available.

Accessories


	B 1.0	14/ 1 1 1 / 1
Unit name	Description	Weight (g)
U-25LBD	LBD filter slider	20
U-25IF550	IF550 filter slider	20
U-25ND6	ND filter	20
U-25ND25	ND filter	20
U-25FR	Frost filter slider	20
U-25L42	UV-cut filter	20
U-PO3	Polarizer slider for reflected light	71
U-POTP3	Polarizer slider for reflected light	71
	with tint plate	
U-AN360-3	360° rotatable analyzer slider	79
U-AN	Analyzer slider for reflected light	50
U-DICR	DIC slider for reflected light	130
U-DICRH	DIC slider for reflected light	130
	(high resolution type)	
U-DICRHC	DIC slider for reflected light	130
	(high contrast type)	





* Combine SZX-TLGAD when using fiber illumination.

Weight: 3.1kg




Reflected light illuminator for BF


U-KMAS

Very compact reflected light illuminator with reduced depth.

Accessories

Unit name	Description	Weight (g)
U-25LBD	LBD filter slider	20
U-25IF550	IF550 filter slider	20
U-25ND6	ND filter	20
U-25ND25	ND filter	20
U-25FR	Frost filter slider	20
U-25L42	UV-cut filter	20
U-PO3	Polarizer slider for reflected light	71
U-POTP3	Polarizer slider for reflected light	71
	with tint plate	
U-AN360-3	360° rotatable analyzer slider	79
U-AN	Analyzer slider for reflected light	50
U-DICR	DIC slider for reflected light	130
U-DICRH	DIC slider for reflected light	130
	(high resolution type)	
U-DICRHC	DIC slider for reflected light	130
	(high contrast type)	





Weight: 1.2kg

Lamp housings

Various different lamp housings are available, for use with different light sources: choose to suit the intended purpose.


U-LH75XEAPO 75W xenon apo lamp housing

(148.5) 30.2 65 65 65 65 65 65 65 65 66 793 88 99 169 (depth dimension for installation) 180.5

Cable length 2,000mm Accepted lamp: UXL-75XB Weight: 3.1kg

*Power supply unit (AH2-RX-T or U-RX-T200) and power cable (UYCP) are necessary for 75W xenon lamp housing. These items are sold separately. AH2-RX-T: dimensions 120(W)x290(D)x186(H), weight approx. 4kg/U-RX-T200 (for EU countries): dimensions 115(W)x195(D)x260(H), weight approx. 3kg Note: Supplied by Olympus Life and Material Science Europa GmbH and its business partners.


U-LH100HGAPO 100W mercury apo lamp housing U-LH100HG 100W mercury lamp housing



* Power supply unit (BH2-RFL-T3 or U-RFL-T200) and power cable (UYCP) are necessary for 100W mercury lamp housings. These items are sold separately.

BH2-RFL-T3: dimensions 120(W)x290(D)x225(H), weight approx 5kg/
U-RFL-T200 (for EU countries): dimensions 150(W)x295(D)x200(H), weight approx. 4.8kg

U-LH100-3/U-LH100IR/U-LH100L-3 100W halogen lamp housings



^{*} External power supply (TH4-100 or TH4-200) and power cable (UYCP) are necessary for 100W halogen lamp housings. These items are sold separately. For TH4-100/200 installation dimensions, refer to the next page.


Lamp housing accessories

For the 100W halogen lamp, the external power supply TH4-100/200 with an intensity adjustment switch and an ON/OFF switch, both are located close to the operator's hand, are provided. All Olympus reflected light illuminators can be used with fiber illumination.

TH4-100/200 External power supply

© OLYMPUS TH4-100 | 82 | 14.5 | 200 | Weight: 2.2kg


TH4-HS Hand switch




U-RMT Extension cord




U-RCV DF converter for BX-URA2



U-LGAD
Fiber adapter for reflected light observation




SZX-TLGAD Transmitted light guide adapter




Light guide mount hole ø12 Weight: 135g

LG-PS2* Light source




*The types of model varies by country in use. Weight: 1.6kg

LG-SF Light guide



U-DULHA Double lamp house adapter




OBSERVATION TUBES

Widefield trinocular observation tubes

Trinocular observation tubes with widefield of view. Compatible with F.N. 22.

U-TR30-2/
Widefield binocular tube
U-TR30IR
Widefield binocular tube for IR

U-ETR-4 Widefield erect image trinocular tube



Unit: mm

Name	Field Number (F.N.)	Inclination angle (degree)	Interpupillary distance (mm)	Light path selector (eyepiece/video port)	Observation image	Weight (g)
U-TR30-2	22	30	50-76	100/0, 20/80, 0/100	Inverted	1600
U-TR30IR	22	30	50-76	100/0, 0/100	Inverted	1600
U-ETR-4	22	30	50-76	100/0, 0/100	Erect	1900

^{*}Length marked with an asterisk (*) may vary according to interpupillary distance. The distance for figure shown is 62mm.

Single port tube with lens

When the visual observation is not needed and only video observation is required, a single port tube with a built-in telan lens can be attached directly to the video port.

U-TLU Single port tube with lens U-TLUIR Single port tube with lens for IR



Weight: 350g


[•] For attachable video camera adapters, refer to video camera adapters system diagram page (pages 5-6).

OBSERVATION TUBES


Super widefield trinocular observation tubes

Trinocular observation tubes with super widefield of view. Compatible with F.N. 26.5.


U-SWTR-3 Super widefield trinocular tube



U-SWETR Super widefield erect image trinocular tube



MX-SWETTR Super widefield erect image tilting trinocular tube



Name	Field Number (F.N.)	Inclination angle (degree)	Interpupillary distance (mm)	Light path selector (eyepiece/video port)	Observation image	Weight (g)
U-SWTR-3	26.5	24	50-76	100/0, 20/80, 0/100	Inverted	2300
U-SWETR	26.5	24	50-76	100/0, 0/100	Erect	4200
MX-SWETTR	26.5	0-42	50-76	100/0, 0/100	Erect	4200

^{*}Length marked with an asterisk (*) may vary according to interpupillary distance. The distance for figure shown is 62mm.


INTERMEDIATE TUBES & ACCESSORIES

Intermediate tubes

Various accessories for various observation need.


U-CA Magnification changer

Provides 1x, 1.2x, 1.6x and 2x intermediate magnifications.




U-TRU Trinocular intermediate attachment

Intermediate attachment which divides the light path, allowing attachment of both digital and video cameras.



U-ECA Magnification changer 2x


Provides 1x and 2x intermediate magnifications.




Weight: 1.3kg

INTERMEDIATE TUBES & ACCESSORIES

U-DP Dual port Use this intermediate tube to divide the light path.





Weight: 1kg


Light path selector by mirror unit




Transmitted side port: side port = 100:0


Transmitted side port: side port = 70:30 (with use of U-MBF3)

U-EPA2 Eyepoint adjuster Raises eyepoint by 30mm.






U-DP1xC Dual port 1x Combine with U-DP to obtain a 1x image.



Weight: 500g


U-APT Arrow pointer Projects an arrow into the field of view.



EYEPIECES/FILAR MICROMETER EYEPIECE

Eyepieces


Eyepieces for UIS2 optical system.




Filar micrometer eyepiece

U-OSM

Used for precise measurement in the field of view.





Weight: 580g

Eyepiece	Magnification 10X, erect image (inverted when used with erect image observation tube), F.N. 14. Diopter adjustment range: ±5 1/m. Provided with rubber eye shade.
Measuring scale	Scale lines graduated in increments of 1mm in the entire 10mm length. Shift of scale lines: 1mm per rotation of the shift ring, the circumference of which is divided into 100 graduations.
Measuring range	10mm/objective lens magnification
Compensation limit for objective lens magnification tolerance	±5% by combined use of the zoom compensation ring and the provided stage micrometer. Compensation ring clamping screw. Magnification compensation scale.
Actual size	Actual size (mm) = Measured value (mm) Objective lens magnification
Repeatability	Repeatability error $\pm \frac{0.007}{A}$ mm (A \cdots Objective lens magnification)
Accuracy	*Measuring error (A ··· Objective lens magnification: L ··· Measured length in mm) ±[(0.0002×A+0.002) L + \frac{0.007}{A}] mm
	Unit: mm

REVOLVING NOSEPIECES

U-D6RE

Sextuple revolving nosepiece with slider slot for DIC


Revolving nosepieces for BF objective lenses

Choose from following 6 types. For motorized nosepieces, refer to motorized unit page.


U-5RE-2 Quintuple revolving nosepiece


U-D6RE-ESD
Sextuple revolving nosepiece with slider slot for DIC with ESD treatment

U-D7RE Septuple revolving nosepiece with slider slot for DIC



Weight: 520g






U-P4RE Centerable quadruple revolving nosepiece with slider slot for DIC

(114.4)
76.4

Weight: 1kg

U-P6RE
Centerable sextuple revolving nosepiece with slider slot for DIC



Insert the DIC dummy when not using the DIC slider

REVOLVING NOSEPIECES


Revolving nosepieces for BF/DF objective lenses


Choose from following 3 types. Use of adapter to mount BF objectives (BD-M-AD) enables attachment of brightfield objective lenses. For motorized nosepieces, refer to motorized unit page.


U-5BDRE Quintuple revolving nosepiece for BF/DF

U-D5BDRE Quintuple revolving nosepiece for BF/DF with slider slot for DIC

U-D6BDRE
Sextuple revolving nosepiece for BF/DF
with slider slot for DIC/
U-P5BDRE
Centerable quintuple revolving nosepiece








Insert the DIC dummy when not using the DIC slider

BD-M-AD Adapter to mount BF objectives




VIDEO CAMERA ADAPTERS


C-mount video camera ports

Allows direct attachment of a C mount video camera. Four types are provided: 0.63x, 0.5x, 0.35x and 0.25x. All models feature a focus adjustment function

U-TV0.25xC C-mount video port with 0.25x lens

U-TV0.35xC-2 C-mount video port with 0.35x lens







Weight: 1.2kg


Weight: 100g

U-TV0.5xC-3 C-mount video port with 0.5x lens

U-TV0.63xC C-mount video port with 0.63x lens







				Unit: mm
]Video camera adapter (Projection lens)	Projection magnifications	Projection area (F.N.)		
		2/3" CCD	1/2" CCD	1/3" CCD
U-TV1x-2	1x	11	8	6
U-TV0.63xC	0.63x	17.5	12.7	9.5
U-TV0.5xC-3	0.5x	22	16	12
U-TV0.35xC-2	0.35x	_	22	17.1
U-TV0.25xC	0.25x	_	_	24

Practical field of view (mm) = Projection area (Field Number)
Objective lens magnifications


Projection area

Focus the video camera adapter to prevent defocusing the eyepiece image and defocusing by magnification switching. Generally, the video camera adapter is focused by switching to a low magnification after focusing at a high magnification objective lens.

VIDEO CAMERA ADAPTERS

Video camera mount adapters


Allows attachment to video cameras with C, Bayonet, Sony and F mounts. Use with the U-TV1x-2. Focus by amount of screwing into U-TV1x-2.



Video camera port

This port can be attached directly to the trinocular observation tube as well as to the single port tube with lens.

U-TV1x-2 Video port 1x



Weight: 150g


Motorized units

Various motorized units, perfect for automation of equipment, are available.

BX-RLAA+U-D6REMC+U-LH100-3


Motorized BF/DF reflected light illuminator+motorized Nomarski DIC sextuple revolving nosepiece+100W halogen lamp housing

Enables motorized exchange of objective lenses, selection between brightfield and darkfield observations as well as aperture diaphragm closing/opening. The BX-UCB control unit has an RS232C connector, allowing control via a PC. For method of attaching illuminator, refer to page 24.




BX-RFAA Motorized universal reflected light illuminator

Reflected light fluorescence illuminator with simultaneous attachment of six mirror units. Incorporates motorized mirror unit changeover and shutter.




U-D5BDREMC Motorized quintuple BD revolving nosepiece with slider slot for DIC U-D6REMC Motorized sextuple revolving nosepiece with slider slot for DIC U-P5REMC


Motorized centerable quintuple revolving nosepiece with slider slot for DIC




BX-UCB Control unit

Motorized units including motorized illuminator and auto focus unit can be totally controlled from BX-UCB






U-HSTR2 Hand switch




BX-REMCB Control box for motorized nosepiece and BF/DF illuminator

BX-RLAA and U-D5BDREMC/U-D6REMC/U-P5REMC can be controlled from U-HSTR2, or direct from the computer keyboard via an RS232C connector.

* BX-RFAA and U-D5BDREM/U-D6REM combination not applicable.




U-ACAD4515 AC adapter for BX-REMCB




^{*} Extension cord U-RMT (1700mm) should be used to connect the lamp housing (U-LH100-3) to the BX-UCB.

U-AFA1M Active auto focus unit




U-FWR Motorized reflected filter wheel

Accomplish maximum 6 filter position exchange




BXFMA-F Motorized illumination with power focus

A motorized microscope unit for integration with your equipment. Motorized operations such as revolving nosepiece up/down, objective lens switching, aperture diaphragm open/close, and brightfield/darkfield switching are accomplished with this component. Several microscopic operations are totally controlled from an external unit by combining this component with an auto focus unit.




* Consult your Olympus dealer about the mounting dimensions.

U-FH Focus adjustment knob unit



U-IFFH Focus adjustment knob interface




Weight: 1450g

DEEP ULTRAVIOLET OBSERVATION SYSTEM


Deep ultraviolet observation system

This module adds a deep ultraviolet (248nm) optical system to a general microscope. An ultra-high resolution observation is executed by using an extremely short wavelength ray.


U-UVF248IM UV248 compatible intermediate tube



U-UVF2FB/5FB UV quartz light guide



U-UVF248LB+U-LH80HGXE UV248 compatible light source box + Mercury Xenon lamp housing



1. Field Number (F.N.) and Practical Field of View

The field number (F.N.) is referred to as the diaphragm size of eyepiece in mm unit which defines the image area of specimen. The diaphragm diameter actually seen through eyepiece is known as the practical field of view (F.O.V.) which is determined by the formula:


$$F.O.V. = \frac{Eyepiece F.N.}{Objective lens magnification} (mm)$$

2. Working Distance (W.D.)

The distance between the front edge of the objective lens and the specimen surface (with the surface of the cover glass in case of the cover glass objective lens) when the specimen is focused.

3. Parfocal Distance

It is the distance between the objective lens mounting plane and the specimen. In UIS2/UIS objective lenses, the parfocal distance is designed at 45mm.



For parfocal distance of the LCPLFLN-LCD series objective lenses, refer to the appropriate objective lens page.

4. Relationship between the objective lens's focal length and magnifications

Indicated magnifications of UIS2/UIS objective lenses are the values when the focal length of the tube lens is 180 mm.

$$M_{(ob)} = \frac{Focal length of tube lens}{f}$$

 $M_{\text{(ob)}}$: Objective lens magnification f: Objective lens's focal length

5. Total Magnification

5.1 Observation through eyepiece (binocular observation)

 $M_{\text{(bino)}} = M_{\text{(ob)}} \times M_{\text{(oc)}}$

 $M_{\text{(bino)}}$: Total magnification for binocular observation

 $M_{\text{(ob)}}$: Objective lens magnification $M_{\text{(oc)}}$: Eyepiece magnification

5.2 Video monitor observation

Total magnification for video monitor

 $M_{\text{(video monitor)}} = M_{\text{(ob)}} \times M_{\text{(video camera adapter)}} \times Monitor magnification*$

M(video monitor): Total magnification on the video monitor

 $M_{\text{(ob)}}$: Objective lens magnification

 $M_{\mbox{(video camera adapter)}}$: Projected magnification for video camera

adapter including photo eyepiece

(refer to Figure 1)

Practical field of view for video monitor observation

Practical field of view for video monitor observation = $\frac{\text{Image device size *}}{M_{\text{(ob)}} \times M_{\text{(video camera adapter)}}}$

M_(ob): Objective lens magnification

 $M_{\mbox{(video camera adapter)}}$: Projected magnification for video camera adapter including photo eyepiece

(refer to Figure 1 for projected magnifications)

Figure 1 Video camera adapter and projection magnifications

Video camera adapter (Projection lens)	Projection magnifications		
U-TV1x-1 + video camera mount adapters	1x		
'			
U-TV0.63xC	0.63x		
U-TV0.5xC-3	0.5x		
U-TV0.35xC-2	0.35x		
U-TV0.25xC	0.25x		

Figure 2 Imaging device size

Camera format	Diagonal	Horizontal	Vertical	
1/3"	6.0mm	4.8mm	3.6mm	
1/2"	8.0mm	6.4mm	4.8mm	
2/3"	11.0mm	8.8mm	6.6mm	

The above table is for standard image device sizes Check your device size for precise calculation.

Figure 3 Imaging device size and monitor magnifications

Camera format	Monitor size (diagonal)				
	9"	12"	14"	21"	27"
1/3"	38.1x	50.8x	59.2x	84.6x	114.1x
1/2"	28.6x	38.1x	44.5x	63.5x	85.7x
2/3"	20.8x	27.7x	32.3x	46.2x	62.3x

Example

What is total magnifications for video monitor when objective lens is 50x, video camera adapter U-TV0.5xC and 2/3" video camera are used?

^{*} Refer to Figure 3 for "Monitor magnification"

^{*} Refer to Figure 2 for image device size

•Total magnification on the video monitor:

 $M_{\text{(ob)}}=50\times$, $M_{\text{(video camera adapter)}}$ is 0.5× from Figure 1 and monitor magnification is 46.2× from Figure 3.

 $M_{\text{(monitor observation)}} = M_{\text{(ob)}} \times M_{\text{(video camera adapter)}} \times \text{monitor magnification}$ = $50 \times 0.5 \times 46.2 = 1155 \times$

•Practical filed of view for video observation(horizontal side): $M_{\text{(ob)}}=50\times$, $M_{\text{(video camera adapter)}}$ is 0.5x from Figure 1 and horizontal side of 2/3" imaging device is 8.8mm from Figure 2

$$\begin{array}{ll} \mbox{Practical field of view} \\ \mbox{for video observation} \end{array} = \frac{\mbox{Image device size}}{\mbox{M}_{(ob)} \times \mbox{M}_{(video \mbox{ camera adapter})}} \\ = \frac{\mbox{8.8 (mm)}}{\mbox{50} \times \mbox{0.5}} = 352 \mu m \end{array}$$

6. Numerical Aperture (N.A.)


The numerical aperture is a key factor to the performance of objective lens (resolving power, focal depth and brightness). The N.A. is determined by the following formula:

$$N.A.= n \times sin\theta$$

- n=Refraction rate of the medium between specimen and objective lenses. (Air: n=1, oil: n=1.515)
- θ: Angle which is made by the optical axis and refraction of the light farthest from the center of the lens.

The visual field brightness (B) of the microscope is determined by the following formula in relation to the objective lens magnification (M). The larger the N.A. and the lower the objective magnification, brightness will increase in the factor of the second power.

$$B \propto \frac{N.A.^2}{M^2}$$



7. Resolving Power

The resolving power of an objective lens is measured by its ability to differentiate two lines or points in an object. The greater the resolving power, the smaller the minimum distance between two lines or points that can still be distinguished. The larger the N.A., the higher the resolving power.

Resolving power formula

The following formula is generally used for determing resolution.

$$\epsilon = 0.61 \times \frac{\lambda}{\text{N.A.}} \text{ (Reyleigh formula)}$$

 λ : Wavelength or radiation in use (λ =0.55µm is used for visible light)

N.A.: Objective lens N.A.

Example

MPLFLN100×(N.A.=0.90), λ =0.55 μ m

$$\epsilon = 0.61 \times \frac{\lambda}{N.A.} = \frac{0.3355}{N.A.} = \frac{0.3355}{0.90} = 0.37 \mu m$$

8. Focal depth of Microscope

The focal depth refers to the depth of the specimen layer which is in sharp focus at the same time, even if the distance between the objective lens and the specimen plane is changed when observing and shooting the specimen plane by microscope. As human eyes are individually different in the ability of their focus adjustment, each person's perception of the focal depth varies. At present, the Berek formula is generally used, because it gives a focal depth value that often coincides with that obtained through experiments.

Focal depth formula

Visual observation (Berek formula)

$$\pm$$
 D.O.F.= $\frac{\omega \times 250,000}{N.A.\times M} + \frac{\lambda}{2(N.A.)^2} (\mu m)$

D.O.F.: Depth Of Focus

- w: Resolving power of eyes 0.0014 (when optical angle is 0.5 degrees)
- M: Total magnification (objective lens magnification x eyepiece magnification)

→ ± D.O.F. =
$$\frac{350}{N.A. \times M}$$
 + $\frac{0.275}{N.A.^2}$ (λ=0.55μm)

This indicates that the focal depth becomes smaller as the numerical aperture becomes larger.

Example

With MPLFLN100x(N.A.=0.90), WHN10x:

$$\pm$$
 D.O.F. = $\frac{350}{0.90 \times 1,000} + \frac{0.275}{0.81} = 0.39 + 0.34 = 0.73 \mu m$

Video camera

In the case of a video camera, the focal depth will vary according to number of pixels of CCD, optical magnification, and numerical aperture. The above-mentioned formula is used as a rough guide only.

9. Aberrations

A difference between an ideal image and an actual image that passes through an optical system is called an "aberration."

9.1 Requirements for Ideal Image Formation

The following three requirements must be satisfied to form an image with no aberration, or an ideal image.

- (i) All the light rays coming from a single point and passing through an image formation optical system converge on a single point.
- (ii) Image points, which correspond to object points on the same plane perpendicular to the optical axis, are present on the same plane.
- (iii) The planar shape of an object and the planar shape of an image that are on the same plane perpendicular to the optical axis have a similarity relation.

Figure 9-1 Requirements for Ideal Image Formation


Object Image plane (ii) (iii)

In an actual optical system, however, it is very difficult to strictly meet the requirements for ideal image formation and this causes "aberrations" that interfere with image forming performance.

9.2 Classification of Aberrations

Aberrations that interfere with image forming performance are classified as shown below in Figure 9-2.

Seidel's aberration = "Expansion of a point image" + "Curvature of image plane" + "Deformation"




Types (1) to (3) correspond to "expansion of a point image" that goes against requirement (i) for ideal image formation in Figure 9-1. Type (4) corresponds to "curvature of image plane" that goes against requirement (ii) in Figure 9-1. Type (5) corresponds to "deformation" that goes against requirement (iii) in Figure 9-1. Types (6) and (7) correspond to "color blur" of images caused by

characteristics of glass materials used for the optical system. "Expansion of a point image" can also be expressed by "wavefront aberration" that regards the light as "waves" and takes account of the phase to include the influence of diffraction.

(1) Spherical aberration


When light rays coming out of an axial object point enter a lens, the light rays with a larger numerical aperture (N.A.) are subjected to stronger refraction power and cross the optical axis in positions with larger differences from the ideal image formation position. The aberration caused this way by different image forming positions due to differences in N.A. of axial light rays is called "spherical aberration." ("Spherical aberration" is proportional to the cube of N.A.)



It is said that objective lenses with larger N.A. have better resolution but worsen spherical aberration. Our advanced design and manufacturing techniques have realized good optical performance even with large numerical aperture.

(2) Coma aberration

Even though spherical aberration is compensated to be very small, there are cases where light rays coming out of an off-axis object point are not condensed to a single point on the image plane but generate asymmetric blur just like a comet leaving traces. This is called coma aberration.



(3) Astigmatism

Even though a lens is compensated for spherical aberration and coma aberration, there are cases where an image of an off-axis object point is not focused to a single point but separated to a concentric line image and a radial line image. This is called "astigmatism." When astigmatism is present, a point image blurs vertically and horizontally, before and after the focus position.

Figure 9-5 Astigmatism and Change in Spot Shape in Different Focus Positions

(a) (b) (c)


(4) Field curvature

An image plane of an object on a plane perpendicular to an optical axis does not always become a plane perpendicular to the optical axis, but it generally becomes a curved plane. This symptom is called "field curvature."

When field curvature is present, the image is more displaced as it becomes closer to the periphery of the visual field. Therefore, when the center of an image is brought into focus, blur occurs in the peripheral areas of the image. To bring the entire image, including the periphery, into clear focus, it is necessary to adequately compensate for this type of aberration.

(5) Distortion

When there is no similar relation between a planar shape on an object and a shape on the image plane, this is called "distortion." When distortion is present, a square image appears in a shape of a barrel or pin-cushion as shown in Figure 9-6.



The microscope optical system contains some distortion. When distortion is present, it can bring erroneous results of shape measurements. When a microscope is used for precision measurements, pay close attention to this aberration, for example, by providing it with an aberration compensation function.

(6) Chromatic aberration


Glasses used for optical systems have different refractive indexes depending on the wavelength. This causes differences in focal length between wavelengths and generates displacement of image forming position. This phenomenon is called "chromatic aberration," which is sometimes subdivided into axial displacement on the optical axis, called "axial chromatic aberration" (or lateral chromatic aberration) and displacement on the image plane, called "chromatic

aberration of magnitude."


Many special glass materials are used, e.g., for apochromats (MPlanApo in Olympus), to eliminate chromatic aberration in a wide range from violet light (g-rays with wavelength of 435 nm) to red light (c-rays with wavelength of 656 nm).

9.3 Wavefront Aberration

Since a long time ago, aberrations have been used in "geometric optics," which considers light as "light rays." Microscope optical systems are often used for observation of very small specimens at a wavelength level, and sometimes adopt "wave optics," which regards light as "waves" and handles the phase information, taking account of the influence of diffraction. In such a case, "wavefront aberration" is used for evaluation. As shown below, when requirements for ideal imaging are satisfied in a microscope optical system, the spherical wavefront (spherical waves) coming from a single point on an object (specimen) is converted to plane waves through an ideal objective lens. The plane waves are converted to spherical waves through an ideal tube lens, and condensed to a single point. The wavefront of these waves is called the "ideal wavefront."



Based on the figure indicated for (1) spherical aberration, the behavior of the wavefront in an optical system that has an aberration is described below.



A difference (a degree of disagreement) between the ideal wavefront and the actual wavefront shown above is called "wavefront aberration."



•OLYMPUS CORPORATION has obtained ISO9001/14001.

Specifications are subject to change without any obligation on the part of the manufacturer.

