
# Anleitung RyecoCam 4000 Digi-Mikroskop

R-FHD-4000-1001

**ryeco**  
microscopes  
+ swiss quality

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# Anleitung RyecoCam 4000 Digi-Mikroskop

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Camera Control Panel	Function	Function Description
	<b>Snap</b>	<b>Capture</b> image from the current video window
	<b>Record</b>	<b>Record</b> video from the current video window
	<b>Auto Exposure</b>	When <b>Auto Exposure</b> is checked, the system will automatically adjust exposure time according to the value of exposure compensation
	<b>Target</b>	Available when <b>Auto Exposure</b> is checked. Slide to left or right to adjust <b>Target</b> according to the current video brightness to achieve proper brightness value
	<b>Exposure Time</b>	Available when <b>Auto Exposure</b> is unchecked. Slide to left or right to reduce or increase exposure time, adjusting brightness of the video
	<b>Gain</b>	Adjust <b>Gain</b> to reduce or increase brightness of video. The Noise will be reduced or increased accordingly
	<b>Red</b>	Slide to left or right to decrease or increase the proportion of <b>Red</b> in RGB on video
	<b>Green</b>	<b>Green</b> is base for reference and cannot be adjusted
	<b>Blue</b>	Slide to left or right to decrease or increase the proportion of <b>Blue</b> in RGB on the video
	<b>White Balance</b>	<b>Auto: White Balance</b> adjustment according to the window video; <b>Manual:</b> Slide the <b>Red</b> or <b>Blue</b> to manually set the video <b>White Balance</b> ; <b>ROI:</b> Set the <b>White Balance</b> according to the <b>ROI</b> . The <b>ROI</b> can be resized and moved;
	<b>Sharpness</b>	Adjust <b>Sharpness</b> level of the video window
	<b>Saturation</b>	Adjust <b>Saturation</b> level of the video window
	<b>Gamma</b>	Adjust <b>Gamma</b> level of the video. Slide to the right side to increase gamma and to the left to decrease gamma.
	<b>Contrast</b>	Adjust <b>Contrast</b> level of the video. Slide to the right side to increase contrast and to the left to decrease contrast.
	<b>DC</b>	For <b>DC</b> illumination, there will be no fluctuation in light source so no need for compensating light flickering
<b>AC(50HZ)</b>	Check <b>AC(50HZ)</b> to eliminate flickering "strap" caused by 50Hz illumination	
<b>AC(60HZ)</b>	Check <b>AC(60HZ)</b> to eliminate flickering "strap" caused by 60Hz illumination	
<b>Default</b>	Set all the settings in the <b>Camera Control Panel</b> to default values	

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Icon	Function	Icon	Function
	<b>Zoom In</b> the Video Window		<b>Zoom Out</b> the Video Window
	<b>Reset Zoom</b> to Original 1X		<b>Horizontal Flip</b>
	<b>Vertical Flip</b>		<b>Video Freeze</b>
	<b>Compare Image</b> with the current video		Display <b>Cross Line</b>
	<b>Browse Images and Videos in the SD Card</b>		<b>Settings</b>
	<b>Check the Version of XCamView</b>		

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Icon	Function
	<b>Float/ Fix</b> switch of the <b>Measurement Toolbar</b>
<input checked="" type="checkbox"/> Visible	Define measuring object in <b>Show up/ Hide</b> mode
Nanometer (nm) ▾	Select the desired <b>Measurement Unit</b>
4X ▾	Choose the same <b>Magnification</b> as the microscope to ensure accuracy of measurement result when measurement unit is not in Pixel unite
	<b>Object Select</b>
	<b>Point</b>
	<b>Angle</b>
	<b>Arbitrary Line</b>
	<b>Parallel</b>
	<b>Horizontal Line</b>
	<b>Vertical Line</b>
	<b>Rectangle</b>
	<b>Circle</b>
	<b>Ellipse</b>
	<b>Annulus</b>
	<b>Two Circles and Center Distance</b>
	<b>Arc</b>

## Anleitung RyecoCam 4000 Digi-Mikroskop

	<b>Polygon</b>
	<b>Curve</b>
	<b>Make Calibration to determine the corresponding relation between magnification and resolution, this will establish the corresponding relationship between measurement unit and the sensor pixel size. Calibration needs to be done with the help of a micrometer. For detailed steps of carrying out calibration please refer to ToupView help manual.</b>
	Export the measurement information to <b>CSV file(*.csv)</b>
	<b>Delete All the Measurement Objects</b>
	<b>Setting</b>
	<b>Exit from Current Measurement Mode</b>
	When the measurement ends, left-click on a single measuring object and the <b>Object Location &amp; Properties Control Bar</b> will show up. The icons on the control bar mean <b>Move Left, Move Right, Move Up, Move Down, Color Adjustment</b> and <b>Delete</b> .

When user left-clicks Display/Hide button  on Measurement Toolbar, the Measurement Toolbar will be fixed. In this case the Camera Control Panel will not pop up automatically even if moving mouse cursor to the left side of the video window. Only when users left click the "X" button on the Measurement Toolbar to exit from the measurement mode will they be able to doing other operations in the Camera Control Panel, the Auto Focus Control Panel or the Camera Control Toolbar.

When a specific measuring object selected during the measuring process, the Object Location and Attributes Control Bar  will appear for changing the object location and properties of the selected objects.

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