RyecoCam FHD Kamera Anleitung

FULL HD RyecoCam no. R-FHD-2000-1201

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RyecoCam FHD Kamera Anleitung

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Camera Control Panel	Function	Function Description
	Snap	Capture image from the current video window
	Record	Record video from the current video window
	Auto Exposure	When Auto Exposure is checked, the system will automatically adjust exposure time according to the value of exposure compensation
	Target	Available when Auto Exposure is checked. Slide to left or right to adjust Target according to the current video brightness to achieve proper brightness value
	Exposure Time	Available when Auto Exposure is unchecked. Slide to left or right to reduce or increase exposure time, adjusting brightness of the video
	Gain	Adjust Gain to reduce or increase brightness of video. The Noise will be reduced or increased accordingly
	Red	Slide to left or right to decrease or increase the proportion of Red in RGB on video
	Green	Green is base for reference and cannot be adjusted
	Blue	Slide to left or right to decrease or increase the proportion of Blue in RGB on the video
	White Balance	Auto: White Balance adjustment according to the window video; Manual: Slide the Red or Blue to manually set the video White Balance ; ROI: Set the White Balance according to the ROI . The ROI can be resized and moved;
	Sharpness	Adjust Sharpness level of the video window
	Saturation	Adjust Saturation level of the video window
	Gamma	Adjust Gamma level of the video. Slide to the right side to increase gamma and to the left to decrease gamma.
	Contrast	Adjust Contrast level of the video. Slide to the right side to increase contrast and to the left to decrease contrast.
	DC	For DC illumination, there will be no fluctuation in light source so no need for compensating light flickering
	AC(50HZ)	Check AC(50HZ) to eliminate flickering "strap" caused by 50Hz illumination
AC(60HZ)	Check AC(60HZ) to eliminate flickering "strap" caused by 60Hz illumination	
Default	Set all the settings in the Camera Control Panel to default values	

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Icon	Function	Icon	Function
	Zoom In the Video Window		Zoom Out the Video Window
	Reset Zoom to Original 1X		Horizontal Flip
	Vertical Flip		Video Freeze
	Compare Image with the current video		Display Cross Line
	Browse Images and Videos in the SD Card		Settings
	Check the Version of XCamView		

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Icon	Function
	Float/ Fix switch of the Measurement Toolbar
<input checked="" type="checkbox"/> Visible	Define measuring object in Show up/ Hide mode
Nanometer (nm) ▾	Select the desired Measurement Unit
4X ▾	Choose the same Magnification as the microscope to ensure accuracy of measurement result when measurement unit is not in Pixel unite
	Object Select
	Point
	Angle
	Arbitrary Line
	Parallel
	Horizontal Line
	Vertical Line
	Rectangle
	Circle
	Ellipse
	Annulus
	Two Circles and Center Distance
	Arc

RyecoCam FHD Kamera Anleitung

	Polygon
	Curve
	Make Calibration to determine the corresponding relation between magnification and resolution, this will establish the corresponding relationship between measurement unit and the sensor pixel size. Calibration needs to be done with the help of a micrometer. For detailed steps of carrying out calibration please refer to ToupView help manual.
	Export the measurement information to CSV file (*.csv)
	Delete All the Measurement Objects
	Setting
	Exit from Current Measurement Mode
	When the measurement ends, left-click on a single measuring object and the Object Location & Properties Control Bar will show up. The icons on the control bar mean Move Left, Move Right, Move Up, Move Down, Color Adjustment and Delete .

1) When user left-clicks Display/Hide button  on the Measurement Toolbar, the Measurement Toolbar will be fixed. In this case the Camera Control Panel will not pop up automatically even if moving mouse cursor to the left side of the video window. Only when users left click the "X" button on the Measurement Toolbar to exit from the measurement mode will they be able to doing other operations in the Camera Control Panel, the Auto Focus Control Panel or the Camera Control Toolbar.

2) When a specific measuring object selected during the measuring process, the Object Location & Attributes Control Bar  will appear for changing the object location and properties of the selected objects.

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