
Binocular tube 30°/20 with integrated HD IP camera for Primotech

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MICROSCOPY + METROLOGY SERVICES
Suisse made.

RYF AG
Bettlachstrasse 2 · CH-2540 Grenchen
Tel. +41 32 654 21 00
ryfag@ryfag.ch · www.ryfag.ch

English

Knowledge of this manual is required for the operation of the instrument. Would you therefore please make yourself familiar with the contents of this manual and pay special attention to hints concerning safe operation of the instrument. The specifications are subject to change; the manual is not covered by an update service. Unless expressly authorized, forwarding and duplication of this document, as well as utilization and communication of its contents are not permitted. Violations will entail an obligation to pay compensation. All rights reserved in the event of granting of patents or registration of a utility model.

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Issued by: Carl Zeiss Microscopy GmbH
Carl-Zeiss-Promenade 10
07745 Jena, Germany

microscopy@zeiss.com
www.zeiss.com/microscopy

Number of this manual: 425538-7144-008
Date of issue: Version 11, 07/2021

CONTENTS

	Page
1	Notes on Instrument Safety 4
1.1	General safety notes 4
1.2	Instrument safety and EMC 5
1.3	Unpacking, transportation, storage 5
1.4	Disposal 5
1.5	Use 6
1.6	Notes on warranty 8
1.7	Intended Use 8
2	Description 9
2.1	System overview 9
3	Start-Up and Operation 10
3.1	Controls and interfaces of Primotech with binocular tube 30°/20 with integrated High Definition (HD) Internet Protocol (IP) camera 10
3.2	Converting the microscope 11
3.2.1	Changing the tube 11
3.3	Binocular tube 30°/20 with integrated HD IP camera 12
3.3.1	Control and functional elements of the binocular tube with integrated camera 12
3.3.2	Connecting the Binocular tube 30°/20 with integrated HD IP camera to the stand 14
3.3.3	Ports and their possible uses 15
3.3.4	Configuring the integrated HD IP camera in the binocular tube 19
4	Care and Troubleshooting 23
4.1	Instrument care 23
4.2	Care of the binocular tube 30°/20 with integrated HD IP camera 23
4.3	Troubleshooting 24
5	Appendix 26
5.1	Technical data 26

1 NOTES ON INSTRUMENT SAFETY

1.1 General safety notes

Please read this Operating Manual carefully before starting up the microscope.

If you need supplementary information, contact our Carl Zeiss Service or an authorized agency.

In addition to this manual, please observe the detailed information given in the following operating manuals: safety information Primotech and Instruction Manual Primotech.

To ensure safe operation and troublefree function of the microscope, strictly observe the precautions and warnings given in this manual.

These are marked herein as follows:

**CAUTION**

This symbol indicates a possible hazard to the user of the instrument.

**CAUTION**

Hot surface!

**CAUTION: LED radiation**

LED Risk Group 1 to IEC 62471, LED radiation will be emitted.

Never look into the LED beam of the illuminating device – either with or without optical instruments.

**ATTENTION**

This symbol indicates a possible hazard to the instrument or system.

**ATTENTION**

Disconnect the plug-in power unit from line power before opening the microscope!

**NOTE**

This symbol refers you to advice that you must observe under all circumstances.

1.2 Instrument safety and EMC

The Primotech microscopes have been designed, produced and tested in compliance with the standards EN 61010-1 (IEC 61010-1) "*Safety requirements for electrical equipment for measurement, control and laboratory use*".

The Primotech microscope complies with Directives 2014/35/EU (Low Voltage Directive), 2014/30/EU (EMC), RoHS Directive 2011/65/EU (including Directive 2015/863 and their amending of Annex II to directive 2011/65/EU) and carries the  mark.

Radio interference suppression in compliance with EN 55011 Class A.

Noise immunity in compliance with EN 61326-1.

The instruments have to be disposed of in compliance with the WEEE Directive 2012/19/EU.

1.3 Unpacking, transportation, storage

Please observe the following safety notes for unpacking, transportation and storage of the microscope:

- The microscope is supplied packed to commercial standards in a plastic case with cardboard packaging; use the original packaging only for any transportation.
- Retain the original packaging for longer storage or return to the manufacturer.
- When unpacking the equipment, verify that all parts specified on the delivery note are present.

- Keep transportation and storage temperatures as specified in Technical Data.
- Set up the microscope on a stable worktable with solid and smooth tabletop.
- Do not touch optical surfaces to avoid fingerprints.

Risk of burns due to hot surface on the underside of the microscope during operation and up to 10 minutes after power off.

1.4 Disposal

Please observe the following safety notes for the disposal of the microscope:

Defective microscopes should not be disposed of with household waste; dispose of them in compliance with the provisions of the law.

The manufacturer of the device is under the legal obligation to take back defective devices.

Batteries of the battery supply unit should not be disposed of with household waste; dispose of them in compliance with the provisions of the law.

1.5 Use

The microscopes including their original accessories must not be used for microscopic techniques other than those described in this Operating Manual.

Please observe the following safety notes when using the microscopes:

The manufacturer cannot assume any liability for other applications, including those of individual modules or single components. This also applies to any service or repair work that is not carried out by authorized service personnel. In case of non-compliance, all warranty claims shall be forfeited.

Opening of the device is only allowed to accordingly instructed specialists or to the Service staff.

Do not operate the devices and their accessories included in the delivery in potentially explosive areas, in the presence of volatile anesthetics or combustible solvents, such as alcohol, benzene or similar chemicals.

Dirt and dust may impair the performance of the devices. They must therefore be protected from such influences to the greatest possible extent and covered with the dust cover when not in use. Always check whether the devices are switched off before you cover them (blue power-on light is off).

The device may only be operated by trained personnel who are aware of the possible dangers involved with microscopy and the particular application. The microscope may only be operated on a stable, solid, smooth and flame resistant surface.

The microscope is a high-precision instrument that can be impaired in its performance or even be destroyed when handled improperly.

The microscope is equipped with a plug-in power unit allowing line voltages to be used in the range between 100 and 240 V ± 10 %, 50 / 60 Hz, without any need for changing the voltage setting on the instrument.

The plug-in power unit meets the requirements of protection class II (with protective insulation). If its casing is damaged, put the plug-in power unit out of operation. The microscope may be operated only with the supplied plug-in power unit.

If it is noted that protection measures are no longer effective, the instrument must be switched off and safeguarded against inadvertent operation. Please contact a Zeiss service agency or the Carl Zeiss Microscopy Service to repair the instrument.

- Always disconnect the power cable before opening the instrument and changing the lamp or LED.

- Wait for the lamp to cool down before replacing it and do not touch the new bulb to avoid fingerprints.

- The instrument may be opened by accordingly instructed specialists or service staff only.
- The operation of the instrument in explosion-risk environments is not allowed.

When using immersion oil, read in any case the safety data sheet.

Immersion oil irritates the skin. Avoid any contact with skin, eyes and clothing.

After skin contact, wash the oil off with plenty of water and soap.

After eye contact, immediately rinse the eye with plenty of water for at least five minutes.

If the irritation persists, consult a medical specialist.

Proper disposal of immersion oil: Take care to ensure that immersion oil does not enter surface water or the sewage system.

The microscope is not equipped with special devices for the protection from substances that are corrosive, potentially infectious, toxic, radioactive, or other substances that might be hazardous to health. Make sure to observe all legal regulations, particularly the relevant national accident prevention regulations when handling such substances.

- Before any transportation of the instrument, switch it off and let it cool down (hot surface on the underside of the instrument).

- Operate the device only on a hard, non-combustible base.
- The plug-in power unit must not get in contact with moisture.

Blocked or covered ventilation slots may lead to heat accumulation that may damage the device and, in extreme cases, cause a fire. Always keep the ventilation slots clear and make sure that no objects are pushed into or fall through the ventilation slots.

The microscope can only be safely disconnected from line power, when the plug-in-power supply is disconnected from line power.

1.6 Notes on warranty

The Primotech microscope and its original accessories are to be used for the microscopy procedures described in this operating manual only. The manufacturer cannot assume any liability for other applications.

Please note the following information on warranty for the microscopes

- The manufacturer guarantees that the microscope is free of material or manufacturing defects when delivered.
- Possible defects must be notified to us immediately and steps be taken to minimize damage.
- If notified of such a defect, the manufacturer is obligated to rectify it at its discretion, either by repairing the instrument or delivering an intact replacement.
- No guarantee is provided for defects caused by natural wear (wearing parts in particular) and improper use.
- The instrument manufacturer shall not be liable for damage caused by faulty operation, negligence or any other tampering with the microscope, particularly the removal or replacement of microscope components, or the use of accessories from other manufacturers.

Unauthorized tampering with the instrument shall lead to a forfeit of all warranty claims.

1.7 Intended Use

Primotech microscopes with Binocular tube 30°/20 with integrated HD IP camera are all-purpose light microscopes primarily designed for industrial applications such as:


- Metallurgy
- Electronics
- Geoscience

Primotech microscopes are also designed for use in education.

2 DESCRIPTION

2.1 System overview

English

For detailed information regarding the system overview of microscope stand, please refer to section 2, "System Overview" in Primotech User Manual (430055-9144-008).

3 START-UP AND OPERATION

3.1 Controls and interfaces of Primotech with binocular tube 30°/20 with integrated High Definition (HD) Internet Protocol (IP) camera


Fig. 1 Controls of Primotech with binocular tube 30°/20 with integrated HD IP camera (425538-9020-000)

Legend to Fig. 1:

- 1 Snap button
- 2 Infrared sensor for remote control
- 3 LED status display
- 4 SD card slot
- 5 Menu button
- 6 AWB button (trigger automatic white balance)
- 7 Increase brightness button
- 8 Reduce contrast button
- 9 Reduce brightness button
- 10 Increase contrast button


Fig. 2 Interfaces on reverse side of Primotech

Legend to Fig. 2:

- 1 HDMI port
- 2 12 V plug-in power unit connection
- 3 Fast Ethernet port
- 4 USB 2.0 port

3.2 Converting the microscope

Unplug the plug-in power unit from line power before converting the microscope.


3.2.1 Changing the tube

For detailed information regarding changing the tube see section *Changing components* in the Primotech Manual (430055-9144-008).

3.3 Binocular tube 30°/20 with integrated HD IP camera

3.3.1 Control and functional elements of the binocular tube with integrated camera

- 1 Snap button
- 2 Infrared sensor for remote control
- 3 LED status display
- 4-9 Function buttons, explanation see Table 1

Do not operate the camera in a potentially explosive environment! Operation of the camera in a potentially explosive environment, in the presence of volatile anesthetics or flammable solvents such as alcohol, benzene or similar may result in explosions or deflagrations.

- **Briefly** press the **Snap** button (Fig. 3/1) to trigger an exposure.
- Press the **Snap** button for **longer** to start a video recording.
- During a video recording **briefly** press the **Snap** button to stop recording.
- Press the **Menu** button (Fig. 3/4) to open the OSD menu.

Fig. 3 Control elements of the binocular tube with integrated camera

Pos. No.	Standard function	Function in OSD (On Screen Display) menu
4	Activate menu	Exit menu
5	AWB (trigger automatic white balance)	OK (apply value)
6	Increase brightness	Right
7	Reduce contrast	Down / reduce value
8	Reduce brightness	Left
9	Increase contrast	Up / increase value


Table 1 Key assignment standard function / function in OSD menu

The visualized images of the integrated camera may only be used for training and research. Direct generation of diagnostic results from these images is not proposed.

The camera has a number of ports:

- SD / SDHC card port (Fig. 6/2) for storing image data, see section 3.3.3.1
- USB 2.0 port (Fig. 4/4) for communication and image data transfer, see section 3.3.3.2
- HDMI port (Fig. 4/1) for image data transfer to a monitor, TV or projector, see section 3.3.3.3
- Fast Ethernet port (Fig. 4/3) for communication and image data transfer, see section 3.3.3.4
- Connection socket for power supply 12 V (Fig. 4/2), see section 3.3.2

- 1 HDMI output
- 2 12 V plug-in power unit connection
- 3 Ethernet port
- 4 USB 2.0 port

Fig. 4 Interfaces on reverse side of the binocular tube 30°/20

The following connecting cables are required for communication and image data transfer and use of the ports:

Item	Name	Figure	Remarks
1	USB 2.0 cable (supplied with A / MINI cable in 3 m length)		Connection between binocular tube 30°/20 with integrated camera and PC
2a	HDMI cable (not supplied)		Connection between binocular tube 30°/20 with integrated camera and monitor, TV or projector
2b	HDMI / DVI adapter cable (not supplied)		Connection between binocular tube 30°/20 with integrated camera and monitor / projector
3	Ethernet cable (two pieces)		Connection between binocular tube 30°/20 with integrated HD IP camera and Microscope Network adapter (MNA). A short network cable is provided as part of the standard delivery.

Table 2 Connections

 If the ZEISS Matscope app is used, a commercially available WLAN router will be required.

The various states of the integrated camera are displayed on the binocular tube 30°/20 by LED.

System status	LED status display / operating mode
Start sequence	LED lights up green
Start sequence completed	LED flashes green
Ready for operation	LED lights up green
Storage process	LED flashes green
Error	LED flashed red for approx. 5 seconds

Table 3 LED status displays


Fig. 5 Connecting the 12 V interconnection cable

3.3.2 Connecting the Binocular tube 30°/20 with integrated HD IP camera to the stand

- Insert the plug (Fig. 5/2) of the 12 V interconnection cable into the connection socket labeled "12V INPUT" (Fig. 5/1) located on the rear of the Binocular tube 30°/20 with integrated HD IP camera.
- Insert the plug on the other end of the 12 V interconnection cable into the connection socket labeled "OUTPUT DC 12V" located on the rear of the stand (Fig. 5/3).
- Connect the plug-in power unit (Fig. 5/4) with a power outlet.

3.3.3 Ports and their possible uses

3.3.3.1 Saving data to the SD card via SD card port

An SD card has a max. capacity of 4 GB, an SDHC (HC = high capacity) up to 32 GB. Prior to use on the PC, the SDHC card must be formatted in the format FAT32.

- Insert an SD / SDHC card (Fig. 6/3) into the card slot (Fig. 6/2) of the binocular tube 30°/20.
- Briefly press the **Snap** button (Fig. 6/1) to trigger an exposure.

The images are saved on the SD / SDHC card in JPG format.

- Press the **Snap** button (Fig. 6/1) for longer to start a video recording. Briefly pressing the **Snap** button once again stops the video recording.


Fig. 6 Inserting SD / SDHC card

Videos are recorded in MPEG-4 advanced video coding format (H.264) and saved as a .h264 file. The videos can be played back using standard software, e.g. with the VLC player.

The visualized images of the integrated camera may only be used for training and research. Direct generation of diagnostic results from these images is not proposed.

3.3.3.2 Data transfer to the PC via USB 2.0 port

Basic imaging acquisition is possible with ZEISS software ZEN or ZEN lite.

- Insert the mini USB plug (Fig. 7/1) into the USB 2.0 port (Fig. 7/2) on the reverse side of the binocular tube 30°/20.
- Insert the USB 2.0 plug (Fig. 7/3) into the PC (Fig. 7/4).
- Connect the 12 V interconnection cable for the tube to get supplied with power by the stand, see section 3.3.2 on page 14.


Fig. 7 Data transfer to the PC via USB port

The visualized images of the integrated camera may only be used for training and research. Direct generation of diagnostic results from these images is not proposed.

The encoding function of Primotech with MNA for the magnification and scaling readout is not valid with ZEN. Only when used in conjunction with iOS software Matscope.

3.3.3.3 Visualization of data via HDMI port

The binocular tube with integrated camera can be connected to a monitor, TV or projector for visualization of the live image data.

- ☞ The data at the HDMI output is in DVI format. Further HDMI functions (e.g. audio, communication between monitor and camera) are not supported.
- ☞ To enable the display the high-resolution images in HD format in this quality, it is recommended that only HD-capable devices be connected.
- Connect the monitor / projector to the binocular tube 30°/20 using the HDMI cable. To do this, insert the HDMI plug (Fig. 8/1) into the jack on the reverse side of the binocular tube 30°/20 (Fig. 8/2) and the monitor / projector (Fig. 8/3). Secure the HDMI plug on the jack on the monitor / projector.


Fig. 8 Connecting the binocular tube 30°/20 to a monitor / projector

- ☞ If a monitor / projector is used, an HDMI / DVI adapter cable (Fig. 8/4) or HDMI / DVI adapter will be required for the reproduction of image data.
- ☞ A default resolution of 720p60 has been set.
- For further camera settings using the on-screen display menu (OSD) see section 3.3.4.
- ☞ Live image data can be documented by photo or video recording, see section 3.3.3.1.
- Connect the 12 V interconnection cable for the tube to get supplied with power by the stand, see section 3.3.2 on page 14.
- ⚠ The visualized images of the integrated camera may only be used for training and research. Direct generation of diagnostic results from these images is not proposed.


Fig. 9 Establishing a network connection

3.3.3.4 Integrate the camera into the network via an Ethernet port

The binocular tube with integrated camera has a fast Ethernet port for communication and image data transfer.

- Integrate the camera into the network via the fast Ethernet port.
To do this, insert one plug (Fig. 9/2) of the short Ethernet cable into the jack on the reverse side of the binocular tube 30°/20 (Fig. 9/1).
- Insert the other plug of the Ethernet cable into the jack of the MNA (Fig. 9/3).
- Insert one plug of the other Ethernet cable into the jack on the side of the MNA (Fig. 9/4) and connect the camera into the network.

A WLAN router will be required to use the built-in camera together with the ZEISS Matscope app. Existing WLAN infrastructure may be used or a separate microscope WLAN set up.

The integrated camera will identify itself automatically to the network (Fig. 9/5) on DHCP (Dynamic Host Configuration Protocol) and will be automatically recognized by Matscope, provided the iPad (Fig. 9/6) is in the same network.

- Connect the 12 V interconnection cable for the tube to get supplied with power by the stand, see section 3.3.2 on page 14.

Consult your network administrator for further guidance.

Further information on the Matscope app can be found at www.zeiss.com/matscope.

For an overview of all ZEISS Microscopy apps visit www.zeiss.com/micro-apps.

The visualized images of the integrated camera may only be used for training and research.

Direct generation of diagnostic results from these images is not proposed.

3.3.4 Configuring the integrated HD IP camera in the binocular tube

- Press the **Menu** button on the tube's control panel.
The OSD menu for adjusting the settings will appear.

The following menus are available for configuring the camera:


Fig. 10 Color menu

AWB	The purpose of white balance is to calibrate the camera to the color temperature of the light source used. The AWB function of the integrated camera supports the "Continuous" / "Push button" modes. In "Continuous" mode the white balance is carried out continuously. In certain specimens this may result in incorrect colors. To ensure that the white balance applies to the light source and not the main color of the specimen, the push button mode can be set. Here the white balance is carried out only once if the AWB / OK button is pressed outside the OSD menu.
Color temp.	This creates a white balance bias towards yellow (left) or blue (right). If AWB (see above) is set to push button, the effect of this regulator is only visible after a white balance has been carried out with the AWB button (outside the OSD menu).
Saturation	The camera's images are in sRGB color space with an 8 bit color depth per channel. The color saturation per RGB value can be adjusted. Extreme left: Gray scales Extreme right: Maximum color saturation
Gamma	The gamma correction of the camera can be adjusted. A choice of Gamma 1 and 0.45 is available.


Fig. 11 Exposure menu

AET / AGC	<p>AET (automatic exposure time) and AGC (auto gain control) automatically ensure the correct brightness of the image when the setting is at "Continuous". The settings are "Once" / "Continuous" / "Off".</p> <p>In some situations (e.g. strong contrasts on the specimen) the automatic setting may cause the image to be under- or overexposed. In this case the automatic setting can be deactivated by switching to "Off".</p> <p>A shorter exposure time with higher amplification (= more noise, but smoother live image) or a longer exposure with lower amplification (= less noise, less smooth live image) can be selected in "Off" mode.</p> <p>In "Once" mode a one-time automatic adjustment will be carried out when confirmed with OK, but then deactivated again.</p>
Auto level	<p>Influences the AET / AGC control when "Continuous" mode is set.</p> <p>Extreme left: very dark image Extreme right: very light image</p>
Exposure	<p>The exposure time can be set within the range 20 μs to ~2 s. The exposure time may be set automatically or manually. Manual setting is possible by selecting the AET / AGC control in "Off" mode.</p> <p>Extreme left: 20 μs Extreme right: 2000 ms (2 s)</p>
Gain	<p>The sensor amplification function can be set within the range 1x to 8x. Sensor gain may be set automatically or manually. Manual setting is possible by selecting the AET / AGC control in "Off" mode.</p> <p>Extreme left: no increased amplification (less image noise, less light sensitivity of sensor)</p> <p>Extreme right: maximum amplification (strong image noise, high light sensitivity of sensor)</p>


Fig. 12 Resolution menu

Preset	<p>The following options are available:</p> <p>720p60: HD resolution, 1280x720, 60 Hz</p> <p>1080p30: Full HD resolution, 1920x1080, 30 Hz</p> <p>A default resolution of 720p60 has been set.</p> <p>Please observe the following: Prior to conversion to the full HD resolution (1080p30), check that the monitor supports this option. If 1080p30 has already been activated and the monitor displays no image, please follow the instructions as described in section 4.3 on page 24.</p>
Video Encoder	<p>The live image and video recordings are influenced by this setting. The following settings are possible:</p> <p>FPS:HALF BR:HIGH: Half frame rate, high data rate (8 Mbit / s)</p> <p>FPS:HALF BR:MID: Half frame rate, medium data rate (4 Mbit / s)</p> <p>FPS:HALF BR: LOW: Half frame rate, low data rate (2 Mbit / s)</p> <p>FPS:FULL BR:HIGH: Full frame rate, high data rate (8 Mbit / s)</p> <p>FPS:FULL BR:MID: Full frame rate, high data rate (4 Mbit / s)</p> <p>FPS:FULL BR:LOW: Full frame rate, low data rate (2 Mbit / s)</p> <p>Half frame rate results in improved image quality with the same data rate.</p>
Sharpening	<p>Value range: 0 to 3.</p> <p>Edge flickering may occur in some HDMI monitors with a default value of 3, but this will disappear with lower values.</p>
Snapshot Res (Res = resolution)	<p>The snapshot resolution for exposures to the SD card can be defined here.</p> <p>The following modes can be set:</p> <ul style="list-style-type: none"> ▪ as monitor on the HDMI port: <ul style="list-style-type: none"> - 1920x1080 (16:9 side ratio), when the camera is in 1080p mode (see default settings) - 1280x720 (16:9 side ratio), when the camera is in 720p mode (see default settings) <p>Advantages: Snap lasts < 1 second, small file sizes due to comparatively low resolution</p> ▪ 5 MP: <ul style="list-style-type: none"> - 2560x1920 (4:3 side ratio) – full sensor resolution, 5 megapixels <p>Advantages: full resolution, wider vertical image field</p> ▪ 2560x1440 <ul style="list-style-type: none"> - 2560x1440 (16:9 side ratio) <p>Advantages: high resolution (3.7 megapixels), same image field as in live image</p>
Snapshot Delay	<p>This setting is useful when images are to be captured using the Snap button on the tube. Because pressing the button causes the microscope to shake, motion blur would be visible in the image with immediate exposure. This effect is eliminated by snapshot delay.</p> <p>The snapshot delay can be set between 0 and 10 seconds.</p>

English

Fig. 13 Setup User menu

Date / Time	Date and time in year-month-day / hour-minute-second format. Each element can be adjusted using the up / down buttons. The new time can be set by pressing the OK button. Date and time are integrated into the folder and file names on the SD card.
Flip image	The following modes are available for image mirroring: FLIP H+V: horizontal and vertical mirroring FLIP V: only vertical mirroring FLIP H: only horizontal mirroring OFF: no mirroring Default setting: the image position on the monitor or iPad corresponds to that in the eyepiece.
Load Settings	The previously saved user settings (USER 1-4) can be loaded in this menu. Selecting "FACTORY DEFAULT" reinstates the default camera settings. After confirmation with OK the camera will be restarted. This may take several seconds. If the OSD menu is not available, the factory settings can be reinstated using the Snap + AWB / OK buttons (keep depressed for two seconds).
Save Settings	The current settings can be saved in a user storage location (USER 1-4). They can be automatically loaded in the "Load Settings" (see above) or "Startup Set" menus (see below) at a later point in time or each time the camera is switched on Please observe the following: If the settings are saved, "Startup Set" is automatically set to this storage location.
Startup Set	This setting determines which user storage location (USER 1-4) is to be loaded each time the camera is switched on.


Fig. 14 Technical Info menu

Version	Current firmware version
MAC	MAC address of the built-in camera
IP	Currently used IP address
Destination IP	Currently used multicast target address
Netmask	Currently used netmask

4 CARE AND TROUBLESHOOTING

4.1 Instrument care

Care of the microscope is restricted to the following operations:

- Cover the instrument with the dust cover after every use.
- Do not install the microscope in a humid room; maximum humidity < 75 %.
- Cover open tubes with the dust caps.
- Remove dust and loose dirt from visible optical surfaces with a brush, blower brush, cotton bud, optics cleaning tissue, or a cotton cloth.
- Remove water-soluble dirt (coffee, cola, etc.) by blowing on it and subsequent wiping it off with a dust-free cotton cloth or a cloth moistened with water to which you may also add a mild detergent.
- Wipe off stubborn oily or fatty dirt (immersion oil, fingerprints) with a cotton bud or a dust-free cotton cloth and the optics cleaning solution L.

The cleaning solution consists of 90 vol% gasoline and 10 vol% isopropanol (IPA). Its individual constituents are also known as:

Gasoline: Surgical spirit, petroleum ether

Isopropanol: 2-propanol,
dimethyl carbinol,
2-hydroxypropane

Clean optical surfaces by polishing in circles starting in the middle and moving to the edges using slight pressure only.

Disconnect the plug-in power unit from the supply line before cleaning it. Make sure that no moisture gets into the plug-in power unit.

For the use in warm and humid climatic zones, all optical components of the microscope are provided with protection against fungus attack.

4.2 Care of the binocular tube 30°/20 with integrated HD IP camera

- The filter glass may only be cleaned with a suitable agent.
- Remove dust from the filter glass with an optical brush, blower brush, cotton swab, optical paper or lint-free cotton cloth.

4.3 Troubleshooting

Problem-solving notes for the use of the built-in camera in conjunction with the network / WLAN and Matscope app are to be found separately at the Zeiss website www.zeiss.com/Matscope.

Problem	Cause	Troubleshooting
LED does not light	The camera does not draw power.	Connect the camera to a plugged-in power supply. Switch on the microscope stand.
LED flashes red	The camera indicates an error.	SD card is full: Insert a new card or delete images from the current card, SD card is write-protected: Push the locking slide on the card to "unlocked". SD card is defective or not formatted: Format SD card; replace SD card if error message does not disappear
Buttons do not react, LED flashes red when a button is pressed.	The buttons are temporarily locked, as the camera is being accessed exclusively from the PC or network.	Control the camera via PC software or Zeiss app, or close the PC software or Zeiss app.
Firmware update not functioning.	For an update an SD card must be inserted and the firmware update must be saved to a specified subfolder on the SD card.	Insert a formatted, unsecured SD card with at least 10 MB free memory space. Be sure to carefully follow the instructions enclosed with the firmware update.
The live image of the camera is not displayed in ZEN.	The camera is not recognized by the Zeiss imaging software.	Several cameras have been connected to the PC: select the correct camera "AxioCamER" from the menu list.
Newly captured images are not visible in Microsoft Windows Explorer on the memory card.	Microsoft Windows Explorer has not updated the contents of the card.	Remove the memory card from the camera, wait for three seconds and re-insert it. Microsoft Windows Explorer will read in the contents of the card again.
The camera "forgets" the date / time.	The buffer battery is empty.	Contact service to replace the buffer battery. Replacement requires special tools and can thus not be carried out by the customer.
The image has severe noise.	The amplification (gain) is set too high.	Reduce the gain and increase the exposure time instead.
The image is too dark or too light.	AET / AGC (automatic exposure time / automatic gain control) has not been activated.	Activate AET / AGC and select "Continuous" mode or manually adjust the exposure time so that the setting is correct for the current light situation.

5 APPENDIX**5.1 Technical data****Dimensions (width x depth x height)**

Stand with binocular tube 30°/20 with integrated HD IP camera for Primotech	approx. 190 mm x 415 mm x 395 mm
---	----------------------------------

Weight

Primotech with binocular tube 30°/20	approx. 8.5 kg
--------------------------------------	----------------

Ambient conditions

Transportation (in packaging): Permissible ambient temperature	-40 to +70 °C
Storage: Permissible ambient temperature Permissible air humidity (no condensation)	+10 to +40 °C max. 75 % at 35 °C
Operation: Permissible ambient temperature Permissible air humidity (no condensation) Atmospheric pressure Operating altitude Degree of pollution	+10 to +40 °C max. 75 % at 35 °C 800 hPa to 1060 hPa max. 2000 m 2

Operating data

Protection class	II
Protection type	IP20
Electrical safety	in compliance with EN 61010-1 (IEC 61010-1) including CSA and UL directives
Pollution degree	2
Overvoltage category	II
Radio interference suppression	in accordance with EN 61326-1
Power supply for Primotech and Binocular tube 30°/20 with integrated HD IP camera:	
Input:	100-240 V AC, 50-60 Hz, max. 0,7 A
Output:	12 V DC, max. 2,5 A
LED illumination	LED risk group 1 acc. to EN 62471:2009 adjustable from 1,5 V to 6 V DC
Binocular tube 30°/20 with integrated HD IP camera for Primotech	
Maximum visual field number (eyepiece)	20
Acquired visual field of the camera	14.25 mm x 10.7 mm (17.82 mm diagonal)
Eyepiece (interpupillary) distance	adjustable from 48 to 75 mm
Viewing angle	30°
Eyepoint height	from 385 mm to 420 mm
Fixed graduation	50 % vis / 50% doc
Optical adaptation	0.4x
Illuminating mirror	with plane surface and spherical surface with f' = 75 mm

HD-CMOS-camera**Sensor-specific data**

Sensor	Micron MT9P031
Sensor size	1/2.5", 5.7 mm x 4.28 mm (7.1 mm diagonal)
Pixel size	2.2 µm x 2.2 µm
Sensor type	1/2.5" CMOS, color
Readout mode	Progressive scan
No. of pixels (H x V), full field	2560 x 1920 pixels active, 5 megapixels
Live image, movie	1920 x 1080 pixels, 30 Fps (H264 max. 16 Mbits / s)
Spectral sensitivity (without IR filter)	400 nm to 700 nm

Signal processing / interface-specific data

Digitization / color depth	24 bit, 3 x 8 bit / pixel
Amplification	0-18 dB
Interfaces	USB 2.0 mini USB plug LAN via RJ 45 jack, 100 Mbit, 12 V plug-in power unit connection, SD card (secure digital) 1-32 GB, slot for SD and SDHC HDMI (1080p/30 or 720p/60)
Remote release	IR sensor
Duo LED	Ready (green), exposure (flashing green), not ready (red), error (flashing red)
Key matrix	White balance, snap, contrast, brightness, menu
Exposure time	10 µs to 2 s

General

Voltage supply connections	via 12 V plug-in power unit through microscope stand
Operating temperature	+5 °C to +45 °C
Relative humidity	max. 80 %, non-condensing

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RYF AG Zweigniederlassung Zürich