

Nikon SMZ 745/745T/445 and SMZ 800N/1270/1270i



Dual Hi-Power Spots



Stand Adapters



LED 40 Ring Lights

Ring light inner \varnothing 66 mm



Working distance 55 -120 mm

P-PS32 stands, column \varnothing 32 mm



619-20-049
Adapter for
P-PS32 stands



619-30-041
619-30-042
619-20-043
85, 200, 300mm High Power Spot
Goose Neck direct mount to base

Objective Adapters

598-32-043

Adapter 66 to 60 mm Obj.
diameter
SMZ 800N/1270/1270i



598-32-050

Adapter with thread
M55x0.75 for
SMZ 745/745T



598-32-049

Adapter with thread
M48x0.75 for SMZ 445




619-30-059


Adapter 66 to 80 mm
Objective diameter
SMZ 18 / 25



P-PS / C-DS stands



619-20-048
Adapter for
P-PS / C-DS stands



598-30-039
Connect fibre optic light source
via Intensilight – light guide to the
Ti diascope illumination pillar

