

Low Cost Basic LED Ring Illuminator CL-16.1



Basic Low Cost LED Ring Illuminator Optika / Ryeco CL16.1:

Basic Ring light Led illuminator, with continuous and external brightness control, external intensity light control, number of LEDs: 144 with total of 500 lumen, four illumination sectors controlled by individual buttons, diameter of the included fixing ring: 60mm, suitable for many stereo microscopes at a reasonable price.

Einfaches LED Ringlicht Optika / Ryeco CL16.1 Basic Ringlicht CL-16.1 LED regelbar Spezifikation: Externe Regelung der Lichtintensität, 144 LEDs, vier Sektoren können einzeln an und aus geschaltet werden, Innendurchmesser der Fixierung: 60mm, passt an viele Stereo Mikroskope, Farbtemperatur (K) 6300.

Éclairage (simple) annulaire LED Optika/ Ryeco CL16.1: Éclairage annulaire LED CL-16.1 pour stéréo microscopes, composé de trois rangées de diodes (144 pièces à 6'300K et 500 lumen) réparties sur la circonférence de l'anneau, permet d'obtenir une lumière blanche homogène. Éclairage simple et réglable avec touches de contrôle, différents modes d'éclairage: circonférence complètement allumée ou seulement un quart ou la moitié. S'adapte (Ø externe à 60mm) sur les sorties objectifs de plusieurs marques grâce aux 3 vis de serrage placées à 120° chacune.


Low Cost Basic LED Ring Illuminator CL-16.1

Key Features Basic CL 16.1:

- Adjustable LED Ring Illuminator with 6300K
- 360° shadow-free illumination
- Solid Machined aluminum housing
- Supplied with 3 tip mounting screws and plastic adapter
- Supplied with separate controller with power supply and cable
- Standard model fits directly all 66mm and 60mm lenses (Adapter included)

Number of LEDs: 144
Inside Ring Diameter: 73 mm
Outside Diameter: 108 mm
Height: 20 mm

Input: 100-240V - 50/60 Hz, 0.5A
Output: 10V – 12V maximum



Ryeco SlimLED
(CA)

Basic Led Ring
CL 16.1 (China)

