



Olympus Adapters - Overview



2 Stand Adapters
 interface to Dual Hi Power Spot
 and Combi Light Guide



2 Lens Adapters


interface to 66 mm inner diameter of
 Photonic ring lights



Fits also to
 35 and 25 mm
 pillars

619-20-049

Adapter for
 pillar 32mm
 SZX2-ST stands
 (SZX 10 SZX 16)



619-20-048

Adapter for
 SZ2 stands
 (SZ 51 / SZ 61)

619-20-056
619-20-033

Back light
 insertable D = 90 / 100 mm



598-32-049

Adapter w thread M48x0.75
 for SZ 51 / SZ 61



598-32-043

Adapter 66 to 60 mm
 SZX 10 / SZX 16



F3000
LED



F3000
LED