


# Adapters - Overview

Only 4 Adapters needed for complete system



**2 Stand Adapters**  
interface to Dual Hi Power Spot  
and Combi Lightguide



Dual Hi  
Power  
Spot



Fits also to  
35 and 25 mm  
pillars

Adapter for  
pillar 32mm  
SZX2-ST stands  
(SZX 10 SZX 16)

**619-20-049**



Adapter for  
SZ2 stands  
(SZ 51 / SZ 61)



**619-20-048**



Combi  
Lightguide



F3000  
LED



**2 Lens Adapters**

interface to 66 mm inner diameter of  
Photonic ringlights



Adapter w. thread M48x0.75  
for SZ 51 / SZ 61

**598-32-049**



Adapter 66 to 60 mm  
SZX 10 / SZX 16

**598-32-043**



F3000  
LED

