



Versatility –
 giving you the freedom to work
 with your favorite equipment.



JENOPTIK ProgRes[®] & PROGRES GRYPHAX[®] microscope cameras

hard- & software compatibility – overview



JENOPTIK microscope cameras

hard- & software compatibility - overview

		PROGRES GRYPHAX [®] software			CapturePro [®] software	
		WIN	MAC	LINUX	WIN	MAC
USB 3.0	PROGRES GRYPHAX [®] SUBRA	X	X	X		
	PROGRES GRYPHAX [®] ARKTUR	X	X	X		
	PROGRES GRYPHAX [®] NAOS	X	X	X		
	PROGRES GRYPHAX [®] RIGEL*	X	X	X		
	PROGRES GRYPHAX [®] KAPELLA*	X	X	X		
	PROGRES GRYPHAX [®] PROKYON*	X	X	X		
USB 2.0	ProgRes [®] CT3 ^{USB}				X	

	ProgRes [®] CT5 ^{USB}				X	
	ProgRes [®] SpeedXTcore 3				X	
	ProgRes [®] SpeedXTcore 5	X			X	
	ProgRes [®] MF ^{USB}				X	
	ProgRes [®] CF ^{USB}				X	
FireWire 400	ProgRes [®] CT3				X	X
	ProgRes [®] C3				X	X
	ProgRes [®] C5				X	X
	ProgRes [®] C7				X	X
	ProgRes [®] MF				X	X
	ProgRes [®] CF				X	X
	ProgRes [®] CF ^{cool}	X	X	X	X	X
	ProgRes [®] MF ^{cool}	X	X	X	X	X
	ProgRes [®] CF _{scan}				X	X
	ProgRes [®] MF _{scan}				X	X
	ProgRes [®] C14 ^{Plus}	X	X	X	X	X

*camera model coming soon

If you need any further information please send an email to:

progres@jenoptik.com

Thank you for your attention!
Your PROGRES GRYPHAX[®] Team

PROGRES
GRYPHAX[®] SERIES

JENOPTIK | Healthcare & Industry
JENOPTIK Optical Systems GmbH
Göschwitzer Straße 25
07745 Jena, Germany