



**fitted with**

**work equipment**  
**connected load**  
**power consumption**  
**glare-free**  
**system of protection**  
**class of protection**  
**usage**  
**luminaire body**

**lamp cover**  
**tubular sections**

**balance of articulation joints**  
**weight (net)**  
**mains lead**

**fastening**

24 x Light-emitting diode  
 light colour neutral white, 4700-5300 K, Color Rendering Index (CRI)=80  
 Energy efficiency category A+  
 electronic transformer  
 100-240 V; 50/60 Hz  
 approx. 31 W  
 conical prismatic screen  
 IP 20  
 II  
 multi-function switch  
 sectional aluminium, anodised, aluminium coloured  
 anodised  
 acrylic clear  
 tubular aluminium section, anodized, aluminium coloured  
 anodised  
 spring  
 approx. 1.9 kg  
 approx. 5 m; mains plug  
 CEE 7/XVI  
 screw-down flange, table clamp (accessory), wall bracket (accessory)



betlachstrasse 2  
 ch-2540 grenchen/suisse  
 Tel. 032 654 21 00  
 fax 032 654 21 09  
 www.ryfag.ch  
**imaging microscopes**




CEE 7/XVI

### illuminance

measuring conditions: d=50cm

Em: 1137 lx  
 Emax: 3053 lx  
 Emin: 159 lx



### luminous intensity

