



MORE LIGHT



System requirements – JENOPTIK GRYPHAX[®] software & USB 3.0 microscope cameras

Beside of the latest software version on the memory card (included in delivery of the JENOPTIK GRYPHAX[®] cameras), we offer a free download after registration.

JENOPTIK GRYPHAX[®] microscope camera software

DOWNLOAD

Note 1: For software and driver installation temporary administrator permission is necessary. Please contact your local IT department in advance!

Note 2: For software operation it is required permission of „full control“ for system folder & subfolders of OS:

Windows OS "C:\ProgramData\Jenoptik\..."
Mac OS "MacHD/Library/Application Support/Jenoptik\..."
Linux OS "/var/lib/Jenoptik/..."

RECOMMENDED SYSTEM REQUIREMENTS (valid as of GRYPHAX version 2.2.0)

Personal Computer	Intel i7 Quad-Core (min. 3.0 GHz) processor 8 GB RAM (dual-channel*)
Operating System	Windows 10 version 1903 or newer – 64 bit MacOS Catalina / Big Sur – 64 bit Linux Ubuntu 18.04 LTS – 64 bit
Data Interface	USB 3.0 integrated over PCI Express V2.0 (Renesas chip set)
Graphic Interface	Dedicated Graphic card equipped with on-board video memory
Monitor Resolutions	1920x1200 for AVIOR & SUBRA & KAPELLA & RIGEL & WEGA & POLARIS & PROKYON 3840x2160 for BETRIA & ALTAIR & ARKTUR & NAOS

* dual-channel memory architecture mandatory to reach full camera performance

MINIMUM SYSTEM REQUIREMENTS

Personal Computer	Intel i5 (Dual-Core) processor or comparable 4 GB RAM (dual-channel*)
Operating System	Windows 10 version 1903 or newer – 64 bit MacOS Catalina / Big Sur – 64 bit Linux Ubuntu 18.04 LTS – 64 bit
Data Interface	USB 3.0 over PCI Express V1.1 (Renesas chip set)
Monitor Resolutions	1280x720 or higher



Wir machen Qualität sichtbar
Nous rendons la qualité visible
Making quality visible

ryfag.ch



MORE LIGHT

NOT SUPPORTED:

JENOPTIK GRYPHAX® cameras on USB2.0 interface! Connection issues & bad images can occur!

JENOPTIK GRYPHAX® cameras under Windows XP or older and 32 bit operating systems!

Learn more about JENOPTIK GRYPHAX® software and the easy workflows and tools.

Watch our video tutorials.

[VIDEO TUTORIALS](#)

HELPFUL? We appreciate your feedback.

[FEEDBACK](#)

We looking forward to hearing from you and remain sincerely,

Your JENOPTIK GRYPHAX® Team

MICROSCOPY • METROLOGY SERVICES
Suisse made



We show quality
Nous rendons la qualité visible
Wir machen Qualität sichtbar



Bettlachstrasse 2 | 2540 Grenchen | Switzerland
phone +41 32 654 21 00 | +41 22 776 82 28
ryfag@ryfag.ch | www.ryfag.ch