

## Universal Ryeco reflected light stand for stereo microscopes

Ryeco RSZ-1090

**ryeco**  
microscopes  
+ swiss quality



Ryeco reflected light stand RSZ-1090 with the Ryf NKL-12 LED ring illuminator

## Universal Ryeco reflected light stand for stereo microscopes

**Ryeco universal incident light stand with a very good price-performance ratio:**

Low-cost stand for various stereomicroscopes such as Ryeco RSZ-45(T), Nikon SMZ 445 / SMZ 745(T), Leica S9, Zeiss 305 / 508 and Olympus SZ51 / SZ61(T). This simple industrial-grade incident light stand with focussing drive has no built-in illumination, but can be supplemented with the most common LED ring light illuminators such as the Ryf RL4b, NKL-12 LED or the most cost-effective variant, the CL-16.1.

Depending on the model, an RGL-4 or an adapter is also required.

The black/white plate, made of Plexiglas, has a diameter of Ø95mm.

The tripod with integrated Z-axis drive is made of aluminium and is painted white (origin Asia).



**Dimensions: W=215mm x D=270mm x H=325mm**

**Ryeco Swiss Warranty & Service Package & Final Assembly**

SAP Order number: RSZ-1090

[www.ryfag.ch](http://www.ryfag.ch)