



User Guide for PROGRES GRYPHAX® Micromanager driver

The “[Micromanager driver](#)” for PROGRES GRYPHAX® USB 3.0 cameras enables user to operate with Micromanager - the open source microscope software solution. The driver delivers essential functionalities of PROGRES GRYPHAX® USB 3.0 cameras. Driver is supported under Windows operating systems.

The GRYPHAX Micromanager driver is part of PROGRES GRYPHAX® software installation (as of version 2.0.0 or newer).

General Preparation:

The Micromanager software version **1.4.23_20170731 – 32 bit** has to be installed before start installation of GRYPHAX Micromanager driver. Recommended USB 3.0 interface card (with Renesas chipset) has to be installed in advance!

Overview:

GRYPHAX camera operating under Micromanager software:

