




Stability –
giving you a reliable research tool
you can count on.



User Guide for TWAIN / DirectX interface for GRYPHAX® USB 3.0 cameras

The “**TWAIN & DirectX driver**” for PROGRES GRYPHAX® USB 3.0 cameras enables user to operate with TWAIN and DirectX supported 3rd party software solutions. The driver delivers basic functionalities of PROGRES GRYPHAX® cameras. Driver is supported under Windows operating systems.

General description:

The GRYPHAX TWAIN & DirectX driver is part of PROGRES GRYPHAX® software installation (as of version 1.1.10 or newer). The driver will be automatically installed and can be used without separate configuration out of the box running a TWAIN or DirectX interface supporting application.

The driver consists of two separate windows. The *preview window* to display the live stream getting from camera and the *settings window*, which contains all essential camera parameters.

The *settings window* is permanently visible, the *preview window* can be deactivated.

